# Go molecular!

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### **Literatur:**

1. Matricardi PM et al. EAACI Molecular Allergy User's Guide Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology. 2016;27 Suppl 23:1-250.
2. Kleine-Tebbe J and Jakob T Editors: Molecular Allergy Diagnostics. Innovation for a Better Patient Management. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6
3. Canonica GW et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics. World Allergy Organ J. 2013 Oct 3;6(1):17.
4. Wickman M. When allergies complicate allergies. Allergy 2005;60(79):14-18.
5. van Hage M et al. ImmunoCAP assays: Pros and cons in allergology. J Allergy Clin Immunol 2017;140:974-7.

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**Literatur:**

1. Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
2. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.
3. Canonica GW et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics. *World Allergy Organ J*. 2013 Oct 3;6(1):17.
4. Sastre J. Molecular diagnosis in allergy. *Clin Exp Allergy* 2010;40(10):1442-1460.
5. Treudler R. and Simon JC. Overview of component resolved diagnostics. *Curr Allergy Asthma Rep* 2013;13(1):110-117.

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**Literatur:**

1. Matricardi PM et al. EAACI Molecular Allergy User's Guide. Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology. 2016;27 Suppl 23:1-250.
2. Kleine-Tebbe J and Jakob T Editors: Molecular Allergy Diagnostics. Innovation for a Better Patient Management. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.
3. Canonica GW et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics. World Allergy Organ J. 2013 Oct 3;6(1):17.
4. Sastre J. Molecular diagnosis in allergy. Clin Exp Allergy 2010;40(10):1442-1460.
5. Treudler R. and Simon JC. Overview of component resolved diagnostics. Curr Allergy Asthma Rep 2013;13(1):110-117. [www.allergen.org](http://www.allergen.org) und [www.allergome.org](http://www.allergome.org)

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### Literatur:

- Nicolaou N et al. Allergy or tolerance in children sensitized to peanut: prevalence and differentiation using component-resolved diagnosis. *J Allergy Clin Immunol* 2010;125:191-197.
- Sicherer SH et al. US prevalence of self-reported peanut, tree nut and sesame allergy: 11 year follow up. *J Allergy Clin Immunol* 2010;125:1322-1326.
- Rona RJ et al. The prevalence of food allergy: a meta-analysis. *J Allergy Clin Immunol* 2007;120(3):638-646.
- Lange L et al. Benefits and limitations of molecular diagnostics in peanut allergy. *Allergo J Int* 2014; 23:158-63.
- Mortz CG et al. The prevalence of peanut sensitization and the association to pollen sensitization in a cohort of unselected adolescents – The Odense Adolescence Cohort Study on Atopic Diseases and Dermatitis (TOACS). *Paediatr Allergy Immunol* 2005;16:501-506.
- Eller E and Bindslev-Jensen C. Clinical value of component-resolved diagnostics in peanut-allergic patients. *Allergy* 2013;68(2):190-194.
- Dang TD et al. Increasing the accuracy of peanut allergy diagnosis by using Ara h 2. *J Allergy Clin Immunol* 2012;129(4):1056-1063.
- Nicolaou N et al. Quantification of specific IgE to whole peanut extract and peanut components in predication of peanut allergy. *J Allergy Clin Immunol* 2011;127(3):684-685.
- Kukkonen AK et al. Ara h 2 and Ara 6 are the best predictors of severe peanut allergy: a double-blind placebo-controlled study. *Allergy* 2015 Oct;70(10):1239-45.
- Rajput S et al. Allergy testing in predicting outcome of open food challenge to peanut. *J Allergy Clin Immunol*. 2018 Jan;141(1):457-458.
- Van Erp FC et al. The IgE and basophil responses to Ara h 2 and Ara h 6 are good predictors of peanut allergy in children. *J Allergy Clin Immunol* 2016. Available online August 8, 2016.
- Klemans RJ et al. The diagnostic accuracy of specific IgE to Ara h 6 in adults is as good as Ara h 2. *Allergy*. 2014 Aug;69(8):1112-4.
- Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
- Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

1. L'Hocine L and Boye J. Allergenicity of soybean: new developments in identification of allergenic proteins, cross-reactivities and hypoallergenization technologies. Crit Rev Food Sci Nutr 2007;47:127-143.
2. Ballmer-Weber B et al. Soy allergy in perspective. Curr Opin Allergy Clin Immunol 2008;8:270-275.
3. Holzhauser T et al. Soybean (Glycinemax) allergy in Europe: Gly m 5 (beta-conglycinin) and Gly m 6 (glycinin) are potential diagnostic markers for severe allergic reactions to soy. J Allergy Clin Immunol 2009;123(2):452-458.
4. Ito T et al. IgE to Gly m 5 and Gly m 6 is associated with severe allergic reactions to soybean in Japanese children. J Allergy Clin Immunol 2010;125;2 Suppl 1:AB88.
5. Kattan DJ and Sampson HJ. Clinical reactivity to soy is best identified by component testing to Gly m 8. J Allergy Clin Immunol Pract. 2015;3(6):970-972.
6. Klemans RJ et al. Components in soy allergy diagnostics: Gly m 2S albumin has the best diagnostic value in adults. Allergy 2013;68:1396-1402.
7. Ebisawa M et al. Gly m 2S albumin is a major allergen with a high diagnostic value in soybean-allergic children. J Allergy Clin Immunol 2013;132:976-978 e1-5.
8. Matricardi PM et al. EAACI Molecular Allergology User's Guide. Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology. 2016;27 Suppl 23:1-250
9. Kosma P et al. Severe reactions after the intake of soy drink in birch pollen-allergic children sensitised to Gly m 4. Acta Paediatr. 2011;100(2):305-306.
10. Kleine-Tebbe J and Jakob T Editors: Molecular Allergy Diagnostics. Innovation for a Better Patient Management. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

- Beyer K et al. Identification of an 11S globulin as a major hazelnut food allergen in hazelnut-induced systemic reactions. *J Allergy Clin Immunol* 2002;110(3):517-523.
- Ortolani C et al. Comparison of results of skin prick tests (with fresh foods and commercial food extracts) and RAST in 100 patients with oral allergy syndrome. *J Allergy Clin Immunol* 1989; 83(3):683-690.
- Faber M et al. Cor a 14: Missing Link in the Molecular Diagnosis of Hazelnut Allergy? *Int Arch Allergy Immunol* 2014;164:200-206.
- Kattan DJ et al. Clinical reactivity to hazelnut may be better identified by component testing than traditional testing methods. *J Allergy Clin Immunol Pract* 2014;2(5):633-634.
- Carraro S et al. COR a 14-specific IgE predicts symptomatic hazelnut allergy in children. *Pediatric Allergy and Immunol* 2016;27(3):322-4.
- Eller E et al. Cor a 14 is the superior serological marker for hazelnut allergy in children, independent of concomitant peanut allergy. *Allergy* 2016;71:556-562.
- Beyer K et al. Predictive values of component-specific IgE for the outcome of peanut and hazelnut food challenges in children. *Allergy* 2015;70:90-98.
- Masthoff L et al. Sensitization to Cor a 9 and Cor a 14 is highly specific for a hazelnut allergy with objective symptoms in Dutch children and adults. *J Allergy Clin Immunol* 2013 Aug;132(2):393-9.
- Brandström J et al. Basophil allergen threshold sensitivity and component-resolved diagnostics improve hazelnut allergy diagnosis. *Clin Exp Allergy* 2015;45(9):1412-8.
- Flinterman AE et al. Lipid transfer protein-linked hazelnut allergy in children from a non-Mediterranean birch-endemic area. *J Allergy Clin Immunol* 2008;121(2):423-428.
- Hansen KS et al. Roasted hazelnuts – allergenic activity evaluated by double-blind, placebo-controlled food challenge. *Allergy* 2003;58(2):132-138.
- Anhoj C et al. Diagnostic evaluation of grass- and birch-allergic patients with oral allergy syndrome. *Allergy* 2001;56(6):548-552.
- Kalyoncu AF et al. Birch pollen related food hypersensitivity: as a para-occupational syndrome. *Allergol Immunopathol (Madr)* 1995;23(2):94-95.
- Bindslev-Jensen C et al. Oral allergy syndrome: the effect of astemizole. *Allergy* 1991;46(8):610-613 .
- Matricardi PM et al. EAACI Molecular Allergy User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
- Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

1. Teuber SS et al. Allergenicity of gourmet nut oils processed by different methods. *J Allergy Clin Immunol* 1997;99(4):502-507.
2. Mew R et al. A retrospect study into the utility of allergen components in walnut allergy. *Ped Allergy and Immunol* 2016;27(7):750-752.
3. Costa J et al. Walnut allergens: molecular characterization, detection and clinical relevance. *Clinical & Experimental Allergy*, 2014(44)319-341.
4. Sato S et al. Jug r 1 sensitization is important in walnut-allergic children and youth. *J Allergy Clin Immunol Pract*. 2017;5(6):1784-1786.
5. Pastorello EA et al. Lipid transfer protein and vicilin are important walnut allergens in patients not allergic to pollen. *J Allergy Clin Immunol* 2004;114(4):908-914.
6. Egger M et al. The role of lipid transfer proteins in allergic diseases. *Curr Allergy Asthma Rep* 2010;20:326-335.
7. Romano A et al. Lipid transfer proteins: The most frequent sensitizer in Italian subjects with food-dependent exercise-induced anaphylaxis. *Clin Exp Allergy* 2012;42(11):1643-1653.
8. Teuber SS et al. Cloning and sequencing of a gene encoding a 2S albumin seed storage protein precursor from English walnut (*Juglans regia*), a major food allergen. *J Allergy Clin Immunol* 1998;101:807-14.
9. Andorf S et al. Association of Clinical Reactivity with Sensitization to Allergen Components in Multifood-Allergic Children. *J Allergy Clin Immunol*. 2017;5(5):1325-1334.
10. Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
11. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

1. Robotham JM et al. Ana o 3, an important cashew nut (*Anacardium occidentale* L.) allergen of the 2S albumin family. *J Allergy Clin Immunol* 2005;115(6):1284-1290.
2. Lange L et al. Ana o 3-specific IgE is a good predictor for clinically relevant cashew allergy in children. *Allergy* 2017;72(4):598-603.
3. Savvatanos S et al. Sensitization to cashew nut 2S albumin, Ana o 3, is highly predictive of cashew and pistachio allergy in Greek children. *J Allergy Clin Immunol* 2015;136(1):192-4.
4. Van der Valk JMP et al. sIgE Ana o 1, 2 and 3 accurately distinguish tolerant from allergic children sensitized to cashew nuts. *Clin Exp Allergy* 2016;47:113-120.
5. Fernandez C et al. Allergy to pistachio: cross reactivity between pistachio nut and other Anacardiaceae. *Clin Exp Allergy* 1995;(25):1254-1259.
6. Parra FM et al. Pistachio nut hypersensitivity: identification of pistachio nut allergens. *Clin Exp Allergy* 1993;23:996-1001.
7. Roux K et al. Tree nut allergens. *Int Arch Allergy Immunology* 2003;131:234-244.
8. Pastorello E et al. Sensitization to the major allergen of Brazil nut is correlated with the clinical expression of allergy. *J Allergy Clin Immunol* 1998;102(6):1021-1027.
9. Andorf S et al. Association of Clinical Reactivity with Sensitization to Allergen Components in Multifood-Allergic Children. *J Allergy Clin Immunol*. 2017;5(5):1325-1334.
10. Brandstrom J et al. IgE to novel citrus seed allergens among cashew-allergic children. *Pediatric Allergy and Immunology* 27 (2016)539-553.
11. O'Sullivan MD and Somerville C. Co-sensitization to orange seed and cashew nut. *Ann Allergy Asthma Immunol* 2011;107:282-3.
12. Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
13. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

1. Rayes H et al. Specific IgE to recombinant protein (Ber e 1) for the diagnosis of Brazil nut allergy. *Clinical & Experimental Allergy* 2015;46, 654-656.
2. Pastorello EA et al. Sensitization to the major allergen of Brazil nut is correlated with the clinical expression of allergy. *J Allergy Clin Immunol* 1998;102:1021-1027.
3. de Leon MP et al. IgE cross-reactivity between the major peanut allergen Ara h 2 and tree nut allergens. *Mol Immunol* 2007;44:463-71.
4. Ewan PW. Clinical study of peanut and tree nut allergy in 62 consecutive patients: new features and associations. *BMJ* 1996;312:1074-8.
5. Pumphrey RSH and Stanworth SJ. The clinical spectrum of anaphylaxis in northwest England. *Clin Exp Allergy* 1996;26:1364-70.
6. Matricardi PM et al. EAACI Molecular Allergy User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
7. Moreno FJ et al. Stability of the major allergen Brazil nut 2S albumin (Ber e 1) to physiologically relevant in vitro gastrointestinal digestion. *FEBS J* 2005;272(2):341-52.
8. Murtagh GJ et al. In vitro stability and immunoreactivity of the native and recombinant plant food 2S albumins Ber e 1 and SFA-8. *Clin Exp Allergy* 2003;33(8):1147-52.
9. Moreno FJ et al. Thermostability and in vitro digestibility of a purified major allergen 2S albumin (Ses i 1) from white sesame seeds (*Sesamum indicum* L). *Biochim Biophys Acta* 2005 Sep 25;1752(2):142-53.
10. Murtagh GJ et al. Stability of recombinant 2 S albumin allergens in vitro. *Biochem Soc Trans* 2001;30(6):913-5.
11. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

1. Brough HA, Caubet JC, Mazon A, Haddad D, Bergmann MM, Wassenberg J et al. Defining challenge-proven coexistent nut and sesame seed allergy: A prospective multicenter European study. *J Allergy Clin Immunol.* 2020;145(4):1231-39.
2. WHO/IUIS (2019). Allergen Nomenclature "Sesamum indicum - All Allergen", WHO/IUIS Allergen Nomenclature Sub-Committee.
3. Pastorello EA, Varin E, Farioli L, Pravettoni V, Ortolani C, Trambaioli C et al. The major allergen of sesame seeds (*Sesamum indicum*) is a 2S albumin. *J Chromatogr B Biomed Sci Appl.* 2001;756(1-2):85-93.
4. Maruyama N, Nakagawa T, Ito K, Cabanos C, Borres MP, Movérare R et al. Measurement of specific IgE antibodies to Ses i 1 improves the diagnosis of sesame allergy. *Clin Exp Allergy.* 2016;46(1):163-71.
5. Tuano KT, Dillard KH, Guffey D, Davis CM. Development of sesame tolerance and cosensitization of sesame allergy with peanut and tree nut allergy in children. *Ann Allergy Asthma Immunol.* 2016;117(6):708-10.
6. Borres M, Maruyama N, Sato S and Ebisawa M. Recent advances in component resolved diagnosis in food allergy. *Allergol Int.* 2016 Oct;65(4):378-387.
7. Sato S, Yanagida N, Ebisawa M. How to diagnose food allergy. *Curr Opin Allergy Clin Immunol.* 2018;18(3):214-21.
8. Foong RX, Dantzer JA, Wood RA, Santos AF. Improving Diagnostic Accuracy in Food Allergy. *J Allergy Clin Immunol Pract.* 2021;9(1):71-80.
9. Yanagida N, Ejiri Y, Takeishi D, Sato S, Maruyama N, Takahashi K et al. Ses i 1-specific IgE and sesame oral food challenge results. *J Allergy Clin Immunol Pract.* 2019;7(6):2084-86.
10. Saf S, Sifers TM, Baker MG, Warren CM, Knight C, Bakhl K et al. Diagnosis of Sesame Allergy: Analysis of Current Practice and Exploration of Sesame Component Ses i 1. *J Allergy Clin Immunol Pract.* 2020;8(5):1681-88.
11. Goldberg MR, Appel MY, Nachshon L, Holmqvist M, Epstein-Rigbi N, Levy MB et al. Combinatorial advantage of Ses i 1-specific IgE and Basophil Activation for diagnosis of Sesame Food Allergy. *Pediatr Allergy Immunol.* 2021 May 5. doi:10.1111/pai.13533. Online ahead of print.
12. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management.* Springer International Publishing Switzerland 2017:152.

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### Literatur:

1. Matricardi PM et al. EAACI Molecular Allergology User's Guide. Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology. 2016;27 Suppl 23:1-250.
2. Faber A et al. IgE-reactivity profiles to nonspecific lipid transfer proteins in a north western European country. J Allergy Clin Immunol. 2017;139(2):679-682.
3. Mothes-Luksch N et al. Pru p 3, a marker allergen for lipid transfer protein sensitization also in Central Europe. Allergy 2017;72:1415-1418.
4. Sastre J. Molecular diagnosis in allergy. Clinical and Experimental Allergy. 2010;40:1442-1460.
5. Scala E et al. Lipid transfer protein sensitization: reactivity profiles and clinical risk assessment in an Italian cohort. Allergy. 2015;70(8):933-43.
6. Pastorello EA et al. Pru p 3-sensitized Italian peach-allergic patients are less likely to develop severe symptoms when also presenting IgE antibodies to Pru p 1 and Pru p 4. Int Arch Allergy Immunol 2011;156:362-372.
7. Klingebiel C. et al. Pru p 7 sensitization is a predominant cause of severe, cypress pollen-associated peach allergy. Clin Exp Allergy 2019;49(4):526-536.
8. Ehrenberg AE et al. Characterization of a 7 kDa pollen allergen belonging to the gibberellin-regulated protein family from three Cupressaceae species. Clin Exp Allergy. 2020 Aug;50(8):964-972 <https://doi.org/10.1111/cea.13675>.
9. Inomata N et al. High prevalence of sensitization to gibberellin-regulated protein (peamaclein) in fruit allergies with negative immunoglobulin E reactivity to Bet v 1 homologs and profilin: Clinical pattern, causative fruits and cofactor effect of gibberellin-regulated protein allergy. J Dermatol. 2017;44(7):735-741.
10. Inomata N et al. Identification of gibberellin-regulated protein as a new allergen in orange allergy. Clin Exp Allergy. 2018;48(11):1509-1520.
11. Tuppo L et al. Pomegranate Cultivars: Identification of the New IgE-Binding Protein Pommaclein and Analysis of Antioxidant Variability. J Agric Food Chem. 2017;65(13): 2702-2710.
12. Kleine-Tebbe J and Jakob T Editors: Molecular Allergy Diagnostics. Innovation for a Better Patient Management. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.
13. Tuppo L et al. Peamaclein – a new peach allergenic protein: similarities, differences and misleading features compared to Pru p 3. Clin Exp Allergy. 2013;43(1):128-140.

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### Literatur:

- Czaja-Bulsa et al. What Do We Know Now about IgE-Mediated Wheat Allergy in Children? *Nutrients* 2017;9(1):35.
- Matricardi PM et al. Primary versus secondary immunoglobulin E sensitization to soy and wheat in the Multi-Centre Allergy Study cohort. *Clin Exp Allergy*. 2008;38(3):493-500.
- Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
- Jones SM et al. Immunologic cross-reactivity among cereal grains and grasses in children with food hypersensitivity. *J Allergy Clin Immunol*. 1995 Sep;96(3):341-51.
- Tatham AS, Shewry PR. Allergens to wheat and related cereals. *Clin Exp Allergy* 2008;38:1712-1726.
- Tanabe S et al. A major wheat allergen has a Gln-Gln-Gln-Pro motif identified as an IgE-binding epitope. *Biochem Biophys Res Commun* 1996;219(2):290-293.
- Battais F et al. Food allergy to wheat: identification of immunoglobulin E and immunoglobulin G-binding proteins with sequential extracts and purified proteins from wheat flour. *Clin Exp Allergy* 2003;33(7):962-970.
- Park HJ et al. Diagnostic Value of the Serum-Specific IgE Ratio of omega-5 Gliadin to Wheat in Adult Patients with Wheat-Induced Anaphylaxis. *Int Arch Allergy Immunol*. 2012;157(2):147-50.
- Morita E et al. Fast omega gliadin is a major allergen in wheat-dependent exercise-induced anaphylaxis. *J Dermatol Sci* 2003;33(2):99-104.
- Palosuo K et al. Rye gamma-70 and gamma-35 secalins and barley gamma-3 hordein cross-react with omega-5 gliadin, a major allergen in wheat-dependent, exercise-induced anaphylaxis. *Clin Exp Allergy* 2001;31(3):466-473.
- Matsuo H et al. Identification of the IgE-binding epitope in omega-5 gliadin, a major allergen in wheat-dependent exercise-induced anaphylaxis. *J Biol Chem* 2004;279(13):12135-12140.
- Palosuo K et al. A novel wheat gliadin as a cause of exercise-induced anaphylaxis. *J Allergy Clin Immunol* 1999;103(5,1):912-917.
- Hofmann SC et al. IgE detection to  $\alpha$ G detection and its clinical relevance in wheat-dependant exercise-induced anaphylaxis. *Allergy* 2012;67:1457-1460.
- Palosuo K et al. Wheat omega-5 gliadin is a major allergen in children with immediate allergy to ingested wheat. *J Allergy Clin Immunol* 2001;108(4):634-638.
- Makela et al. Wheat allergy in children - new tools for diagnostics. *Clin Exp Allergy* 2014;44:1420-1430.
- Ito K et al. IgE antibodies to  $\omega$ -5 gliadin associate with immediate symptoms on oral wheat challenge in Japanese children. *Allergy* 2008;63:1536-1542.
- Nilsson N et al. Wheat allergy in children evaluated with challenge and IgE antibodies to wheat components. *Pediatr Allergy Immunol* 2015;26:119-125.
- Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

- Chokshi NY et al. Molecular diagnosis of egg allergy: an update. *Expert Rev Mol Diagn.* 2015;15(7):895-906.
- Aabin B et al. Identification of IgE-binding egg white proteins: comparison of results obtained by different methods. *Int Arch Allergy Immunol* 1996;109(1):50-57.
- Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology.* 2016;27 Suppl 23:1-250.
- Ando H et al. Utility of ovomucoid-specific IgE concentrations in predicting symptomatic egg allergy. *J Allergy Clin Immunol* 2008;122:583-588.
- Lemon-Mulé H et al. Immunological changes in children with egg allergy ingesting extensively heated egg. *J Allergy and Clin Immunol* 2008;122:977-983.
- Urisu A. Allergenic activity of heated and ovomucoid-depleted egg white. *J Allergy Clin Immunol* 1997;100:171-176.
- Benhamou Senouf AH et al. Native and denatured egg white protein IgE tests discriminate hen's egg allergic from egg-tolerant children. *Pediatr Allergy Immunol* 2015;26:12-17.
- Gray CL et al. Egg sensitization, allergy and component patterns in African children with atopic dermatitis. *Pediatr Allergy Immunol* 2016;27:709-15.
- Bernhisel-Broadbent J et al. Allergenicity and antigenicity of chicken egg ovomucoid (Gal d 1) compared with ovalbumin (Gal d 2) in children with egg allergy and in mice. *J Allergy Clin Immunol* 1994;93:1047-1059.
- Jarvinen KM et al. Specificity of IgE antibodies to sequential epitopes of hen's egg ovomucoid as a marker for persistence of egg allergy. *Allergy* 2007;62:758-765.
- Benhamou AH et al. State of the art and new horizons in the diagnosis and management of egg allergy. *Allergy* 2010;65:283-289.
- Gradman J et al. Relationship between specific IgE to egg components and natural history of egg allergy in Danish children. *Pediatr Allergy Immunol.* 2016 Dec;27(8):825-830.
- Dhanapala P et al. Cracking the egg: An insight into egg hypersensitivity. *Mol Immunol.* 2015;66(2):375-83.
- De Silva C et al. Molecular and immunological analysis of hen's egg yolk allergens with a focus on YGP42 (Gal d 6). *Mol Immunol.* 2016;71:152-60.
- Yanagida N et al. Safety and feasibility of heated egg yolk challenge for children with egg allergies. *Pediatr Allergy Immunol.* 2017;28(4):348-354.
- Hemmer W et al. Update on the bird-egg syndrome and genuine poultry meat allergy. *Allergo J Int.* 2016;25:68-75.
- Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management.* Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

1. Matricardi PM et al. EAACI Molecular Allergology User's Guide. Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology. 2016;27 Suppl 23:1-250.
2. Kaiser C et al. Cow's milk-protein allergy: results of skin-prick test with purified milk proteins. Z Ernahrungswiss 1990;29:122-128.
3. Werfel T et al. Milk-responsive atopic dermatitis is associated with a casein-specific lymphocyte response in adolescent and adult patients. J Allergy Clin Immunol 1997;99(1):124-133.
4. Norgaard A et al. Allergenicity of individual cow milk proteins in DBPCFC-positive milk allergic adults. J Allergy Clin Immunol 1996;97:237.
5. Dupont D et al. Food processing increases casein resistance to simulated infant digestion. Mol Nutr Food Res 2010;54(11):1677-1689.
6. Docena G et al. Identification of casein as the major allergenic and antigenic protein of cow's milk. Allergy 1996;51(6):412-416.
7. Shek LP. Humoral and cellular responses to cow milk proteins in patients with milk- induced IgE-mediated and non-IgE mediated disorders. Allergy 2005;60(7):912-919.
8. Lam HY. Cow's milk allergy in adults is rare but severe; both casein and whey proteins are involved. Clin Exp Allergy 2008;38(6):995-1002.
9. Bloom A et al. Effect of heat treatment on milk and egg proteins allergenicity. Pediatric Allergy and Immunology 2015;25:740-746.
10. Nowak-Wegrzyn AK, et al. Tolerance to extensively heating milk in children with cow's milk allergy. J Allergy Clin Immunol 2008;122(2):342-347.
11. Caubet, JC et al. Utility of casein-specific IgE levels in predicting reactivity to baked milk. J Allergy Clin Immunol 2012;131:222-224.
12. Ito K et al. The usefulness of casein-specific IgE and IgG4 antibodies in cow's milk allergic children. Clin Mol Allergy 2012 Jan 2;10(1):1. doi: 10.1186/1476-7961-10-1.
13. Bartuzi Z et al. Contribution of Molecular Allergen Analysis in Diagnosis of Milk Allergy. Curr Allergy Asthma Rep. 2017;17(7):46.
14. Chatchatee P et al. Identification of IgE and IgG binding epitopes on beta- and kappa-casein in cows milk allergic patients. Clin Exp Allergy 2001;31:1256-62.
15. Chatchatee P et al. Identification of IgE- and IgG-binding epitopes on alpha (s1)-casein: differences in patients with persistent and transient cow's milk allergy. J Allergy Clin Immunol 2001;107:379-83.
16. Cerecedo I et al. Mapping of the IgE and IgG4 sequential epitopes of milk allergens with a peptide microarray-based immunoassay. J Allergy Clin Immunol 2008;122:589-594.
17. Caubet JC et al. Natural tolerance development in cow's milk allergic children: IgE and IgG4 epitope binding. Allergy. 2017 Mar 27. doi: 10.1111/all.13167. [Epub ahead of print]
18. Poza-Guedes P. Role of specific IgE to  $\beta$ -lactoglobulin in the gastrointestinal phenotype of cow's milk allergy. Allergy Asthma Clin Immunol. 2016 Feb 23;12:7.
19. Werfel SJ. Clinical reactivity to beef in children allergic to cow's milk. J Allergy Clin Immunol 1997 99(3):293-300.
20. Martelli A, et al. Beef allergy in children with cow's milk allergy; cow's milk allergy in children with beef allergy. Ann Allergy Asthma Immunol 2002;89(6):Suppl1:38-43.
21. Kleine-Tebbe J and Jakob T Editors: Molecular Allergy Diagnostics. Innovation for a Better Patient Management. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

1. Commins SP et al. Delayed anaphylaxis, angioedema, or urticaria after consumption of red meat in patients with IgE antibodies specific for galactose-alpha-1,3-galactose. *J Allergy Clin Immunol.* 2009 Feb;123 (2):426-33.
2. Saleh H et al. Anaphylactic Reactions to Oligosaccharides in Red Meat: a Syndrome in Evolution. *Clinical and Molecular Allergy* 2012;10(1):5.
3. Kennedy JL et al. Galactose-alpha-1,3-galactose and delayed anaphylaxis, angioedema, and urticaria in children. *Pediatrics* 2013;131:e1545-52.
4. Commins SP et al. Delayed clinical and ex vivo response to mammalian meat in patients with IgE to galactose-alpha-1, 3-galactose. *J Allergy Clin Immunol.* 2014 Jul;134(1):108-15.
5. Caponetto P et al. Gelatin-containing sweets can elicit anaphylaxis in a patient with sensitization to galactose-alpha-1,3-galactose. *J Allergy Clin Immunol Pract.* 2013 May-Jun;1 (3):302-3.
6. Rispens T et al. IgE production to alpha-gal is accompanied by elevated levels of specific IgG1 antibodies and low amounts of IgE to blood group B. *PLoS One.* 2013;8(2):e55566.
7. Morisset M et al. Anaphylaxis to pork kidney is related to IgE antibodies specific for galactosealpha-1,3-galactose. *Allergy.* 2012;67(5):699-704.
8. Fischer J et al. Galactose-alpha-1,3-galactose sensitization is a prerequisite for pork-kidney allergy and cofactor-related mammalian meat anaphylaxis. *J Allergy Clin Immunol* 2014;134:755-759.
9. van Nunen S. Galactose-alpha-1,3-galactose, mammalian meat and anaphylaxis: a worldwide phenomenon? *Curr Treat Options Allergy.* 2014;1(3):262-77.
10. van Nunen SA et al. An association between tick bite reactions and red meat allergy in humans. *Med J Aust.* 2009 May 4;190(9):510-1.
11. Hamsten C et al. Identification of galactose-alpha-1,3-galactose in the gastrointestinal tract of the tick *Ixodes ricinus*; possible with red meat allergy. *Allergy* 2013;68(4):549-52.
12. Hamsten C et al. Red meat allergy in Sweden: association with tick sensitization and B-negative blood groups. *J Allergy Clin Immunol.* 2013 Dec;132(6):1431-4.
13. Villalta D et al. Galactose-alpha-1,3-galactose syndrome: an Italian survey. *Eur Ann Allergy Clin Immunol.* 2017;49(6):263-269.
14. Chung CH et al. Cetuximab-induced anaphylaxis and IgE specific for galactose- $\alpha$ -1,3-galactose. *NEJM* 2008;358 (11):1109-17.
15. Mullins RJ et al. Relationship between red meat allergy and sensitization to gelatin and galactose- $\alpha$ -1,3-galactose. *J Allergy Clin Immunol.* 2012 May;129(5):1334-1342.
16. Matricardi PM et al. EAACI Molecular Allergology User's Guide. Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology. 2016;27 Suppl 23:1-250.
17. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management.* Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

1. Leung NYH, et al. Current immunological and molecular biological perspectives on seafood allergy: A comprehensive review. *Clin Rev Allergy Immunol.* 2014;46(3):180-97.
2. Mariona P et al. Molecular Diagnosis of Shrimp Allergy: Efficiency of Several Allergens to Predict Clinical Reactivity *Journal of Allergy and Clinical Immunol: In Practice* 2015;3(4):521-529.
3. Shanti KN et al. Identification of tropomyosin as the major shrimp allergen and characterization of its IgE-binding epitopes. *J Immunol* 1993;151(10):5354-5363.
4. Yu CJ et al. Proteomics and immunological analysis of a novel shrimp allergen, Pen m 2. *J Immunol* 2003;170:445-453.
5. Ayuso R et al. Myosin light chain is a novel shrimp allergen, Lit v 3. *J Allergy Clin Immunol* 2008;122:795-802.
6. Shiomi K et al. Sarcoplasmic calcium-binding protein: identification as a new allergen of the black tiger shrimp *Penaeus monodon*. *Int Arch Allergy Immunol* 2008;146:91-98.
7. Giuffrida MG et al. Shrimp allergy beyond Tropomyosin in Italy: clinical relevance of Arginine Kinase, Sarcoplasmic calcium binding protein and Hemocyanin. *Eur Ann Allergy Clin Immunol.* 2014 46;(5):172-177.
8. Matricardi PM et al. EAACI Molecular Allergy User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology.* 2016;27 Suppl 23:1-250.
9. Lopata AL and Lehrer SB. New insights into seafood allergy. *Curr Opin Allergy Clin Immunol* 2009;9:270-277.
10. DeWitt AM et al. Recombinant tropomyosin from *Penaeus aztecus* (rPen a 1) for measurement of specific immunoglobulin E antibodies relevant in food allergy to crustaceans and other invertebrates. *Mol Nutr Food Res* 2004;48(5):370-379.
11. Fernandes J. Immunoglobulin E antibody reactivity to the major shrimp allergen, tropomyosin, in unexposed Orthodox Jews. *Clin Exp Allergy* 2003;33:956.
12. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management.* Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

1. Sharp MF et al. Fish allergy: in review. *Clin Rev Allergy Immunol* 2014;46:258-271.
2. Swoboda I, et al. Recombinant carp parvalbumin, the major cross-reactive fish allergen: a tool for diagnosis and therapy of fish allergy. *Allergy* 2002; 57: Suppl 73:79-84.
3. Bugajska-Schretter A et al. Parvalbumin, a cross-reactive fish allergen, contains IgE-binding epitopes sensitive to periodate treatment and Ca<sup>2+</sup> depletion. *J Allergy Clin Immunol* 1998 101:67-74.
4. Lim DL-C et al. Parvalbumin - the major tropical fish allergen. *Pediatr Allergy Immunol* 2008 19:399-407.
5. Bugajska-Schretter A et al. Purification, biochemical, and immunological characterisation of a major food allergen: different immunoglobulin E recognition of the apo- and calcium-bound forms of carp parvalbumin. *Gut* 2000;46(5):661-669.
6. Griesmeier U et al. Expression levels of parvalbumins determine allergenicity of fish species. *Allergy* 2010;65:191-198.
7. Kuehn A et al. Identification of enolases and aldolases as important fish allergens in cod, salmon and tuna: component resolved diagnosis using parvalbumin and the new allergens. *Clin Exp Allergy* 2013;43:811-822.
8. Matricardi PM et al. EAACI Molecular Allergy User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
9. Kuehn A. Fish Allergens at a Glance: Variable Allergenicity of Parvalbumins, the Major Fish Allergens. *Front Immunol*. 2014;5:179.
10. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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## Literatur:

1. Popescu FD. Cross-reactivity between aeroallergens and food allergens. *World J Methodol* 2015 June 26;5(2):31-50.
2. Murray S et al. Preventing severe asthma exacerbations in children: A randomised trial of mite impermeable bedcovers. *Am J Respir Crit Care Med*. 2017 Jul 15;196(2):150-158.
3. Morgan JW et al. Results of Home based Environment Intervention among Urban Children with Asthma. *N Engl J Med* 2004;351:1068-1080.
4. Canonica GW et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics World Allergy Organization Journal 2013;6(1):17. 7.
5. Asero R. Component-resolved diagnosis-assisted prescription of allergen-specific immunotherapy: a practical guide *Eur Ann Allergy Clin Immunol*. 2012;44(5):183-7.
6. Schmid-Grendelmeier P et al. Recombinant allergens – routine diagnostics or still only science? *Der Hautarzt* 2010;61(11):946-953.
7. Focke M et al. Heterogeneity of commercial timothy grass pollen extracts. *Clin Exp Allergy* 2008;38(8):1400-1408.
8. Focke M et al. Molecular composition and biological activity of commercial birch pollen allergen extracts. *Eur J Clin Invest* 2009;39(5):429-436.
9. Brunetto B et al. Characterization and comparison of commercially available mite extracts for in vivo diagnosis. *Allergy* 2010;65(2):184-90.
10. Casset A et al. Varying allergen composition and content affects the in vivo allergenic activity of commercial *Dermatophagoides pteronyssinus* extracts. *Int Arch Allergy Immunol*. 2012;159(3):253-62.
11. Moreno Benitez F et al. Variation in allergen content in sublingual allergen immunotherapy with house dust mites. *Allergy*. 2015;70(11):1413-20.
12. Weghofer M et al. Identification of Der p 23, a peritrophin-like protein, as a new major *Dermatophagoides pteronyssinus* allergen associated with the peritrophic matrix of mite fecal pellets. *J Immunol*. 2013;190(7):3059-67. allergen immunotherapy with house dust mites. *Allergy*. 2015;70(11):1413-20.

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## Literatur:

1. Hilger C et al. Animal Lipocalin allergens. *Curr Allergy Asthma Rep* 2012;12:438-447.
2. Nordlund B et al. IgE antibodies to animal-derived lipocalin, kallikrein and secretoglobulin are markers of bronchial inflammation in severe childhood asthma. *Allergy* 2012;67:661-9.
3. Konradsen JR et al. Allergy to Furry animals: New Insights, diagnostic approaches, and challenges. *J Allergy Clin Immunol* 2015;135:616-25.
4. Werfel SJ. Clinical reactivity to beef in children allergic to cow's milk. *J Allergy Clin Immunol* 1997;99(3):293-300.

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### Literatur:

- Ingram JM et al. Quantitative assessment of exposure to dog (Can f 1) and cat (Fel d 1) allergens: Relation to sensitization and asthma among children living in Los Alamos, New Mexico. *J Allergy Clin Immunol* 1995;96:449-56.
- Bjerg A et al. A population-based study of animal component sensitization, asthma, and rhinitis in schoolchildren. *Pediatr Allergy Immunol.* 2015;26(6):557-63.
- Konradsen JR et al. Severe childhood asthma and allergy to furry animals: Refined assessment and using molecular based allergy diagnostics. *Pediatr Allergy Immunol.* 2014;25:187-192.
- Nordlund B et al. IgE antibodies to animal-derived lipocalin, kallikrein and secretoglobulin are markers of bronchial inflammation in severe childhood asthma. *Allergy* 2012;67:661-669.
- Konradsen JR et al. Allergy to Furry animals: New Insights, diagnostic approaches, and challenges. *J Allergy Clin Immunol* 2015;135:616-25.
- Nagao M et al. Sensitization to secretoglobulin and lipocalins in a group of young children with risk of developing respiratory allergy. *Clin Mol Allergy* 2017;3;15:4. doi: 10.1186/s12948-017-0061-8. eCollection 2017.
- Asarroj A et al. Sensitization to cat and dog allergen molecules in childhood and prediction of symptoms of cat and dog allergy in adolescence: A BAMSE/MeDALL study. *J Allergy Clin Immunol.* 2016;137(3):813-21.
- Drouet M and Sabbah A. The Pork/Cat Syndrome or Crossed Reactivity between Cat Epithelia and Pork Meat. *Monogr Allergy* 1996; 32:164-73.
- Posthumous J et al. Initial description of pork-cat syndrome in the United States. *J Allergy Clin Immunol* 2013;131:924–925 (letter to the editor).
- Hilger C et al. Allergic cross-reactions between cat and pig serum albumin. *Allergy* 1997;52:179-87.
- Canonica GW et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics World Allergy Organization Journal 2013;6(1):17. 7.
- Asero R. Component-resolved diagnosis-assisted prescription of allergen-specific immunotherapy: a practical guide *Eur Ann Allergy Clin Immunol.* 2012;44(5):183-7.
- Schmid-Grendelmeier P et al. Recombinant allergens – routine diagnostics or still only science? *Der Hautarzt* 2010;61(11):946-953.
- Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology.* 2016;27 Suppl 23:1-250.
- Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management.* Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.
- Smith W, O'Neil SE, Hales BJ, et al. Two newly identified cat allergens: the von Ebner gland protein Fel d 7 and the latherin-like protein Fel d 8. *Int Arch Allergy Immunol* 2011;156:159-70.
- Hilger C, Kuehn A, Hentges F. Animal lipocalin allergens. *Curr Allergy Asthma Rep* 2012;12:438-47.
- Apostolovic D et al. The cat lipocalin Fel d 7 and its cross-reactivity with the dog lipocalin Can f 1. *Allergy.* 2016 Oct;71(10):1490-5.

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### Literatur:

- Perzanowski MS et al. Effect of cat and dog ownership on sensitization and development of asthma among pre-teenage children. *Am J Respir Crit Care Med* 2002;166:696-702.
- Nordlund B et al. IgE antibodies to animal-derived lipocalin, kallikrein and secretoglobin are markers of bronchial inflammation in severe childhood asthma. *Allergy* 2012;67:661-9.
- Nicholas C et al. Dog characteristics and allergen levels in the home. *Ann Allergy Asthma Immunol* 2010;105:228-33.
- Konradsen JR et al. Allergy to Furry animals: New Insights, diagnostic approaches, and challenges. *Allergy Clin Immunol* 2015;135:616-25.
- Mattsson L. Molecular and immunological characterization of Can f 4: a dog dander allergen cross-reactive with a 23 kDa odorant-binding protein in cow dander. *Clin Exp Allergy* 2010; 40(8):1276-1287.
- Jakob T et al. Clinical relevance of sensitization to cross-reactive lipocalin Can f 6. *Allergy* 2013;68(5):690-691.
- Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
- Mattsson L et al. Prostatic kallikrein: A new major dog allergen. *J Allergy Clin Immunol* 2009; 123(2):362-368.
- Basagaña M. Involvement of Can f 5 in a Case of Human Seminal Plasma Allergy *Int Arch Allergy Immunol* 2012;159:143-146.
- Kofler L. et al. A case of dog-related human seminal plasma allergy *Eur Ann Allergy Clin Immunol* 2012 Apr;44(2):89-92.
- Schoos AM et al. Precision allergy: Separate allergies to male and female dogs. *J Allergy Clin Immunol Pract*. 2017 Nov-Dec;5(6):1754-1756.
- Canonica GW, et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics *World Allergy Organization Journal* 2013;6(1):17. 7.
- Asero R. Component-resolved diagnosis-assisted prescription of allergen-specific immunotherapy: a practical guide *Eur Ann Allergy Clin Immunol*. 2012;44(5):183-7.
- Schmid-Grendelmeier P et al. Recombinant allergens – routine diagnostics or still only science? *Der Hautarzt* 2010;61(11):946-953
- Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.
- Apostolovic D et al. The cat lipocalin Fel d 7 and its cross-reactivity with the dog lipocalin Can f 1. *Allergy*. 2016 Oct;71(10):1490-5.

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### Literatur:

1. Gawlik et al. Anaphylaxis as a manifestation of horse allergy. WAO Journal 2009;2:185-189.
2. Cosme-Blanco W et al. Anaphylaxis to Horses and Epinephrine Use: Increasing Awareness Among Pediatric Patients and Families. Pediatr Allergy Immunol 2017;28(6):608-610.
3. Roberts G and Lack G. Horse allergy in children. BMJ. 2000;321: 286-287.
4. Konradsen JR et al. Allergy to Furry animals: New Insights, diagnostic approaches, and challenges. Allergy Clin Immunol 2015;135:616-25.
5. Nordlund B et al. IgE antibodies to animal-derived lipocalin, kallikrein and secretoglobulin are markers of bronchial inflammation in severe childhood asthma. Allergy 2012;67:661-9.
6. Matricardi PM et al. EAACI Molecular Allergology User's Guide. Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology. 2016;27 Suppl 23:1-250.
7. Canonica GW et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics World Allergy Organization Journal 2013;6(1):17. 7.
8. Asero R. Component-resolved diagnosis-assisted prescription of allergen-specific immunotherapy: a practical guide Eur Ann Allergy Clin Immunol. 2012;44(5):183-7.
9. Schmid-Grendelmeier P et al. Recombinant allergens – routine diagnostics or still only science? Der Hautarzt 2010;61(11):946-953.
10. Kleine-Tebbe J and Jakob T Editors: Molecular Allergy Diagnostics. Innovation for a Better Patient Management. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

1. Calderon MA et al. House Dust Mite Respiratory Allergy: An Overview of Current Therapeutic Strategies. *J Allergy Clin Immunol Pract.* 2015;3(6):843-55.
2. Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology.* 2016;27 Suppl 23:1-250.
3. Nolte H et al. Major allergen content consistency of SQ house dust mite sublingual immunotherapy tablets and relevance across geographic regions. *Annals of allergy, asthma & immunology: official publication of the American College of Allergy, Asthma, & Immunology.* 2016;117(3):298-303.
4. Weghofer M et al. Identification of Der p 23, a peritrophin-like protein, as a new major Dermatophagoides pteronyssinus allergen associated with the peritrophic matrix of mite fecal pellets. *J Immunol.* 2013;190(7):3059-67.
5. Becker S et al. Real-Life Study for the Diagnosis of House Dust Mite Allergy - The Value of Recombinant Allergen-Based IgE Serology. *Int Arch Allergy Immunol.* 2016;170(2):132-7.
6. Casset A et al. Varying allergen composition and content affects the in vivo allergenic activity of commercial Dermatophagoides pteronyssinus extracts. *Int Arch Allergy Immunol.* 2012;159(3):253-62.
7. Mueller GA et al. Serological, genomic and structural analyses of the major mite allergen Der p 23. *Clin Exp Allergy.* 2016;46(2):365-76.
8. Posa D et al. Evolution and predictive value of IgE responses toward a comprehensive panel of house dust mite allergens during the first 2 decades of life. *J Allergy Clin Immunol* 2017;139:541-549.
9. Resch Y et al. Different IgE recognition of mite allergen components in asthmatic and non asthmatic children. *The Journal of allergy and clinical immunology.* 2015;136(4):1083-91.
10. Canonica GW, et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics *World Allergy Organization Journal* 2013;6(1):17. 7.
11. Asero R. Component-resolved diagnosis-assisted prescription of allergen-specific immunotherapy: a practical guide *Eur Ann Allergy Clin Immunol.* 2012;44(5):183-7.
12. Schmid-Grendelmeier P et al. Recombinant allergens – routine diagnostics or still only science? *Der Hautarzt* 2010;61(11):946-953.
13. Thomas WR. House Dust Mite Allergens: New Discoveries and Relevance to the Allergic Patient. *Current allergy and asthma reports.* 2016;16(9):69.
14. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management.* Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

1. Barber D et al. Understanding patient sensitization profiles in complex pollen areas: a molecular epidemiological study. *Allergy*. 2008 Nov; 63(11):1550-8.
2. Andersson K et al. Characteristics and immunobiology of grass pollen allergens. *International Archives of Allergy & Immunology*. 2003;130(2): 87-107.
3. Hatzler L et al. Molecular spreading and predictive value of preclinical IgE response to Phleum pratense in children with hay fever. *J Allergy Clin Immunol*. 2012 Oct;130(4):894-901 e5.
4. Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
5. Sekerkova A et al. Detection of Phl p 1, Phl p 5, Phl p 7 and Phl p 12 specific IgE antibodies in the sera of children and adult patients allergic to Phleum pollen. *Allergol Int*. 2012 Jun; 61(2):339-46.
6. Tripodi S et al. Molecular profiles of IgE to Phleum pratense in children with grass pollen allergy: Implications for specific immunotherapy. *J Allergy Clin Immunol*. 2012 Mar;129(3): 834-9 e8.
7. Cipriani F et al. Diagnostic relevance of IgE sensitization profiles to eight recombinant Phleum pratense molecules. *Allergy* 2017;Oct 20. doi: 10.1111/all.13338. [Epub ahead of print].
8. Hauser M et al. Panallergens and their impact on the allergic patient. *Allergy Asthma Clin Immunol*. 2010;6(1):1.
9. Schmid-Grendelmeier P. Recombinant allergens – routine diagnostics or still only science? *Der Hautarzt* 2010;61(11):946-953.
10. Focke M et al. (2008) Heterogeneity of commercial timothy grass pollen extracts. *Clin Exp Allergy* 38(8):1400-1408.
11. Walker SM et al. Immunotherapy for allergic rhinitis. *Clin Exp Allergy*. 2011 Sep; 41(9): 1177-200.
12. Valenta R. et al. Component-resolved diagnosis to optimize allergen-specific immunotherapy in the Mediterranean area. *J Investig Allergol Clin Immunol*. 2007;17 Suppl 1:36-40.
13. Canonica GW, et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics World Allergy Organization Journal 2013;6(1):17. 7.
14. Asero R. Component-resolved diagnosis-assisted prescription of allergen-specific immunotherapy: a practical guide *Eur Ann Allergy Clin Immunol*. 2012;44(5):183-7.
15. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

- Hauser M et al. Panallergens and their impact on the allergic patient. *Allergy, Asthma & Clinical Immunology*. 2010;6(1)1.
- Rossi RE et al. Sensitization profiles in polysensitized patients from a restricted geographical area: Further lessons from multiplexed component resolved diagnosis. *Eur Ann Allergy Clin Immunol*. 2011; 43(6):171-175.
- Hauser M et al. Bet v 1-like pollen allergens of multiple Fagales species can sensitize atopic individuals. *Clinical & Exp Allergy*. 2011;41:1804-181.
- Vieths S et al. Current understanding of cross-reactivity of food allergens in pollen. *Ann N.Y Acad Sci*. 2002;964:47-68.
- Schmidt-Andersen MB et al. Identification of European allergy patterns to the allergen families PR-10, LTP and profiling from Rosaceae fruits. *Clin Rev Allerg Immunol*. 2009; 41(1):4-19.
- Valenta R et al. Component-Resolved Diagnosis to Optimize Allergen-Specific Immunotherapy in the Mediterranean area. *J Invest Allergol Clin Immunol*. 2007;Vol 17, supplement 1:88-92.
- Schmid-Grendelmeier P. Recombinant allergens. For routine use or still only science? *Hautarzt*. 2010;61(11):946-53.
- Seckerková A et al. Detection of Bet v 1, Bet v 2 and Bet v 4 specific IgE antibodies in the sera of children and adult patients allergic to birch pollen: evaluation of different IgE reactivity profiles depending on age and local sensitization. *Int Arch Allergy Immunol*. 2011;154:278-285.
- Canonica GW, et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics World Allergy Organization Journal 2013;6(1):17. 7.
- Asero R. Component-resolved diagnosis-assisted prescription of allergen-specific immunotherapy: a practical guide *Eur Ann Allergy Clin Immunol*. 2012;44(5):183-7.
- Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
- Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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## Literatur:

1. García BE et al. Oleaceae-induced pollinosis in an area with exposure to olive and ash trees. *Journal of Investigational Allergology and Clinical Immunology* 2011;21(1):34-37.
2. Castro AJ et al. Pla I 1 and Ole e 1 pollen allergens share common epitopes and similar ultrastructural localization. *J Investig Allergol Clin Immunol* 2007;17 Supplement 1:41-47.
3. Rodríguez R et al. Olive pollen recombinant allergens: value in diagnosis and immunotherapy. *J Investig Allergol Clin Immunol* 2007;17 Suppl 1:4-10.
4. Rossi RE et al. Sensitization profiles in polysensitized patients from a restricted geographical area: Further lessons from multiplexed component resolved diagnosis. *European Annals of Allergy and Clinical Immunology* 2011;43(6):171-175.
5. Barber D et al. Understanding patient sensitization profiles in complex pollen areas: A molecular epidemiological study *Allergy* 2008;63(11):1550-1558.
6. Quiralte J et al. Modelling diseases: The allergens of *Olea europaea* pollen. *J Investig Allergol Clin Immunol* 2007;17 Suppl 1:24-30.
7. Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
8. Hauser M et al. Pan allergens and their impact on the allergic patient. *Allergy, Asthma & Clinical Immunology* 2010; 6(1):1-14.
9. Asturias JA et al. Purification and characterization of Pla a 1, a major allergen from *Platanus acerifolia* pollen. *Allergy* 2002;57(3):221-7.
10. Asturias J et al. The major *Platanus acerifolia* pollen allergen Pla a 1 has sequence homology to invertase inhibitors. *Clin Exp Allergy* 2003;33(7):978-985.
11. Lauer I et al. Identification of a plane pollen lipid transfer protein (Pla a 3) and its immunological relation to the peach lipid-transfer protein, Pru p 3. *Clin Exp Allergy* 2007;37:261-269.
12. Scala E et al. Lipid transfer protein sensitization: reactivity profiles and clinical risk assessment in an Italian cohort. *Allergy* 2015;70:933-943.
13. Arilla MC et al. Quantification of the Major Allergen from Cypress (*Cupressus arizonica*) Pollen, Cup a 1, by Monoclonal Antibody-Based ELISA. *Int Arch Allergy Immunol*. 2004 May;134(1):10-6. Epub 2004 Mar 25
14. Di Felice G et al. Cupressaceae Pollinosis: Identification, Purification and Cloning of Relevant Allergens. *Int Arch Allergy Immunol* 2001;125:280-289.
15. Charpin D et al. Cypress Pollinosis: from Tree to Clinic. *Clin Rev Allergy Immunol*. 2017 Apr 11 doi: 10.1007/s12016-017-8602-y. [Epub ahead of print].
16. Dominguez-Ortega J et al. Prevalence of allergic sensitization to conifer pollen in a high cypress exposure area. *Allergy Rhinol (Providence)*. 2016 Jan 1;7(4):200-206.
17. Caimmi D et al. Epidemiology of cypress pollen allergy in Montpellier. *J Investig Allergol Clin Immunol*. 2012;22(4):280-5.
18. Aceituno E et al. Molecular cloning of major allergen from *Cupressus arizonica* pollen: Cup a 1. *Clin Exp Allergy* 2000;30:1750-1758.
19. WHO/IUIS Allergen Nomenclature Sub-committee. <http://www.allergen.org>.
20. Douladiris N et al. A molecular diagnostic algorithm to guide pollen immunotherapy in Southern Europe: towards component resolved management of allergic diseases. *Int Arch Allergy Immunol* 2013;162:163-172.
21. Asam C et al. Tree pollen allergens - an update from a molecular perspective. *Allergy* 2015;70:1201-1211.
22. Schmid-Grendelmeier P. Recombinant allergens. For routine use or still only science? *Hautarzt*. 2010; 61(11): 946-53.
23. Canonica GW, et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics World Allergy Organization Journal 2013;6(1):17. 7.
24. Asero R. Component-resolved diagnosis-assisted prescription of allergen-specific immunotherapy: a practical guide *Eur Ann Allergy Clin Immunol*. 2012;44(5):183-7.
25. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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## Literatur:

1. Matricardi PM et al. EAACI Molecular Allergy User's Guide. Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology. 2016;27 Suppl 23:1-250.
2. Gadermaier G et al. Allergens of weed pollen: An overview on recombinant and natural molecules. Methods 2014;66:55-66.
3. Hirschwehr R et al., Identification of common allergenic structures in mugwort and ragweed pollen. J Allergy Clin Immunol 1998;101(2 Pt 1):196-206.
4. Asero R et al. Concomitant sensitization to ragweed and mugwort pollen: who is who in clinical allergy? Ann Allergy Asthma Immunol 2014;113:307-313.
5. Fuchs T et al. Natural latex, grass pollen, and weed pollen share IgE epitopes. J Allergy Clin Immunol 1997;100(3):356-64.
6. Helbling A. Food allergy. [German] Ther Umsch 1994;51(1):31-7.
7. Egger M et al. Pollen food syndromes associated with weed pollinosis: an update from the molecular point of view. Allergy 2006;61:461-476.
8. van Toorenenbergen AW et al. Demonstration of spice-specific IgE in patients with suspected food allergies. J Allergy Clin Immunol 1987;79(1):108-13.
9. Jensen-Jarolim E et al. Characterization of allergens in Apiaceae spices: anise, fennel, coriander and cumin. Clin Exp Allergy 1997;27(11):1299-306.
10. Schmid-Grendelmeier P. Recombinant allergens. For routine use or still only science? Hautarzt. 2010; 61(11): 946-53.
11. Canonica GW, et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics World Allergy Organization Journal 2013;6(1):17. 7.
12. Asero R. Component-resolved diagnosis-assisted prescription of allergen-specific immunotherapy: a practical guide Eur Ann Allergy Clin Immunol. 2012;44(5):183-7.
13. Kleine-Tebbe J and Jakob T Editors: Molecular Allergy Diagnostics. Innovation for a Better Patient Management. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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## Literatur:

1. Medrek SK et al. Fungal sensitization is associated with increased risk of life-threatening asthma. J Allergy Clin Immunol Pract. 2007;5:1025-31.
2. Denning DW et al. The link between fungi and severe asthma: a summary of the evidence. Eur Respir J. 2006 Mar;27(3):615-26.
3. Rick EM et al. Allergic Fungal Airway Disease. J Investig Allergol Clin Immunol. 2016;26(6):344-354.
4. Castanhinha S et al. Pediatric severe asthma with fungal sensitization is mediated by steroid-resistant IL-33. J Allergy Clin Immunol. 2015 Aug;136(2):312-22.
5. Moreno A et al. Orthologous Allergens and Diagnostic Utility of Major Allergen Alt a 1. Allergy Asthma Immunol Res. 2016 Sep;8(5):428-37.
6. Gabriel MF et al. *Alternaria alternata* allergens: Markers of exposure, phylogeny and risk of fungi-induced respiratory allergy. Environ Int. 2016 Apr-May;89-90:71-80.
7. Canonica GW et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics. World Allergy Organ J. 2013 3;6(1):17.

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## Literatur:

1. Bush RK and Prochnau JJ. *Alternaria*-induced asthma. *J Allergy Clin Immunol* 2004; 113:227-34.
2. Black PN et al. Sensitivity to fungal allergens is a risk factor for life-threatening asthma. *Allergy*. 2000 May;55(5):501-4.
3. Vailes LD et al. IgE and IgG antibody responses to recombinant Alt a 1 as a marker of sensitization to *Alternaria* in asthma and atopic dermatitis. *Clin Exp Allergy*. 2001 Dec;31(12):1891-5.
4. Downs SH et al. Clinical importance of *alternaria* exposure in children. *Am J Respir Crit Care Med* 2001;164:455-9.
5. Halonen M et al. *Alternaria* as a major allergen for asthma in children raised in a desert environment. *Am J Respir Crit Care Med* 1997;155:1356-61.
6. Halonen M et al. Two subphenotypes of childhood asthma that differ in maternal and paternal influences on asthma risk. *Am J Respir Crit Care Med* 1999;160(2):564-70.
7. Corsico R et al. Prevalence of sensitization to *Alternaria* in allergic patients in Italy *Ann Allergy Asthma Immunol* 1998;80(1):71-6.
8. Dowaisan A et al. Sensitization to aeroallergens among patients with allergic rhinitis in a desert environment. *Ann Allergy Asthma Immunol* 2000;84(4):433-8.
9. Bartra J et al. Sensitization to *Alternaria* in patients with respiratory allergy. *Front Biosci (Landmark Ed)*. 2009 Jan 1;14:3372-9.
10. Kustrzeba-Wójcicka I et al *Alternaria alternata* and its allergens: a comprehensive review. *Clin Rev Allergy Immunol*. 2014 Dec;47(3):354-65.
11. Twaroch TE et al. Carrier-Bound Alt a 1 Peptides without Allergenic Activity for Vaccination Against *Alternaria Alternata* Allergy. *Clin Exp Allergy* 2012; 42(6):966-975.
12. Unger A et al. Clinical testing of recombinant allergens of the mold *Alternaria alternata*. *Int Arch Allergy Immunol* 1999;118:220-1.

13. De Vouge MW et al. Isolation and expression of a cDNA clone encoding an *Alternaria alternata* Alt a 1 subunit. *Int Arch Allergy Immunol* 1996;111:385-95.
14. Postigo I et al. Diagnostic value of Alt a 1, fungal enolase and manganese-dependent superoxide dismutase in the component-resolved diagnosis of allergy to Pleosporaceae. *Clin Exp Allergy*. 2011 Mar;41(3):443-51.
15. Moreno A et al. Orthologous Allergens and Diagnostic Utility of Major Allergen Alt a 1. *Allergy Asthma Immunol Res*. 2016 Sep;8(5):428-37.
16. Gabriel MF et al. *Alternaria alternata* allergens: Markers of exposure, phylogeny and risk of fungi-induced respiratory allergy. *Environ Int*. 2016 Apr-May;89-90:71-80.
17. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

1. Fukutomi Y et al. Serological diagnosis of allergic bronchopulmonary mycosis: Progress and challenges. *Allergology International* 2016;65:30-36.
2. Carsin A et al. *Aspergillus fumigatus* in cystic fibrosis: An update on immune interactions and molecular diagnostics in allergic bronchopulmonary aspergillosis. *Allergy*. 2017;72:1632-1642.
3. Bowyer P et al. Relative reactivity of *Aspergillus* allergens used in serological tests *Medical Mycology* 2006;44:23-28.
4. H. Tanimoto et al. Molecular-based allergy diagnosis of allergic bronchopulmonary aspergillosis in *Aspergillus fumigatus*-sensitized Japanese patients *Clinical & Experimental Allergy* 2015;45, 1790–1800 + erratum *Clin Exp Allergy* 2016;46(2):381.
5. Kurup VP. *Aspergillus* antigens: which are important? *Medical Mycology Supplement* 1 2005;43 (1): s189-S196.
6. Canonica GW, et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics *World Allergy Organ J* 2013;6(1):17. 7.
7. Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
8. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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## Literatur:

1. Spillner E et al. Hymenoptera allergens: from venom to "venome". *Frontiers in immunology* 2014; 5:1-7.
2. Bilo B et al. EAACI Interest Group on Insect Venom Hypersensitivity. Diagnosis of Hymenoptera venom allergy. *Allergy* 2005;60:1339-1349.
3. Bonifazi F et al. and EAACI Interest Group on Insect Venom Hypersensitivity. Prevention and treatment of hymenoptera venom allergy: guidelines for clinical practice. *Allergy* 2005;60:1459-1470.
4. Muller U et al. IgE to recombinant allergens Api m 1, Ves v 1, and Ves v 5 distinguish double sensitization from cross-reaction in venom allergy. *Allergy* 2012;67:1069-1073.
5. Mittermann I et al. Recombinant allergen-based IgE testing to distinguish bee and wasp allergy et al. June 2010. Volume 125, Issue 6, 1300-1307.e3.
6. Michael J et al. Added sensitivity of component-resolved diagnosis in hymenoptera venom-allergic patients with elevated serum tryptase and/or mastocytosis. *Allergy* 2016;71:651-660.
7. Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric Allergy and Immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
8. Kohler et al. Component resolution reveals additional major allergens in patients with honeybee venom allergy. *J Allergy Clin Immunol* 2014;133:1383-9.
9. Korošec P et al. High sensitivity of CAP-FEIA rVes v 5 and rVes v 1 for diagnosis of *Vesputa* venom allergy *J Allergy Clin Immunol*. 2012 May;129(5):1406-8.
10. Grunwald T et al. Molecular cloning and expression in insect cells of honeybee venom allergen acid phosphatase (Api m 3). *Allergy Clin Immunol* 2006;117:848-54.
11. Blank S et al. Api m 10, a genuine *A. mellifera* venom allergen, is clinically relevant but underrepresented in therapeutic extracts. *Allergy* 2011;66:1322-1329.
12. Bonadonna P et al. Clonal mast cell disorders in patients with systemic reactions to Hymenoptera stings and increased serum tryptase levels. *J Allergy Clin Immunol* 2009;123:680-6.
13. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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### Literatur:

- Sutherland MF et al. Specific monoclonal antibodies and human immunoglobulin E show that Hev b 5 is an abundant allergen in high protein powdered latex gloves. *Clin Exp Allergy* 2002;32(4):583-589.
- Rozynek P et al. Cloning, expression and characterization of the major latex allergen prohevein. *Clin Exp Allergy* 1998;28(11):1418-1426.
- Raulf-Heimsoth M et al. Characterization of B- and T-cell responses and HLA-DR4 binding motifs of the latex allergen Hev b 6.01 (prohevein) and its post-transcriptionally formed proteins Hev b 6.02 and Hev b 6.03. *Allergy* 2004;59(7):724-733.
- Vandenplas O et al. The role of allergen components for the diagnosis of latex-induced occupational asthma. *Allergy* 2016;71:840-849.
- Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
- Wagner B, et al. Hev b 7 is a Hevea brasiliensis protein associated with latex allergy in children with spina bifida. *J Allergy Clin Immunol* 2001;108(4):621-627.
- Ebo DG et al. Component-resolved diagnosis from latex allergy by micro-array. *Clin Exp Allergy* 2010;40(2):348-358.
- Ott H et al. Microarrays of recombinant Hevea brasiliensis proteins: A novel tool for the component-resolved diagnosis of natural rubber latex allergy. *J Investig Allergol Clin Immunol* 2010;20(2):129-138.
- Schuler S et al. Microarray-based component-resolved diagnosis of latex allergy: isolated IgE-mediated sensitization to latex profilin Hev b 8 may act as confounder. *Clin Transl Allergy* 2013;3(1):11.
- O'Riordain G et al. Cloning and molecular characterization of the Hevea brasiliensis allergen Hev b 11, a class I chitinase. *Clin Exp Allergy* 2002;32(3):455-62.
- Nettis E Diagnosis of latex allergy: the importance of hev B 11. *Int Arch Allergy Immunol*. 2012;159 (2) 147-8.
- Garnier L et al. Molecular allergens in the diagnosis of latex allergy. *Eur Ann Allergy Clin Immunol* 2012;44(2):73-79.
- Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.

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## Literatur:

1. Jensen-Jarolim E et al. Debates in allergy medicine: Molecular allergy diagnosis with ISAC will replace screenings by skin prick test in the future. *World Allergy Organization Journal*. 2017;10:33.
2. Onell A et al. Allergy testing in children with persistent asthma: comparison of four diagnostic methods. *Allergy* 2017;72: 590-597.
3. Foong R et al. Pilot study: assessing the clinical diagnosis of allergy in atopic children using a microarray assay in addition to skin prick testing and serum specific IgE. *Clin Mol Allergy* 2016; 14:8.
4. Canonica GW, et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics. *World Allergy Organ J* 2013;6(1):17.
5. Matricardi PM et al. EAACI Molecular Allergology User's Guide. *Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology*. 2016;27 Suppl 23:1-250.
6. Kleine-Tebbe J and Jakob T Editors: *Molecular Allergy Diagnostics. Innovation for a Better Patient Management*. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.
7. Patellis A, et al. Multiplex component-based allergen microarray in recent clinical studies *Clinical & Experimental Allergy*. 2016;46:1022-1032.
8. Melioli G et al. Allergenius, an expert system for the interpretation of allergen microarray results. *World Allergy Organ J* 2014;7:15.
9. Panzner P et al. A comprehensive analysis of middle-European molecular sensitization profiles to pollen allergens. *Int Arch Allergy Immunol*. 2014;164:74-82.
10. Huss-Marp J et al. Comparison of molecular and extract-based allergy diagnostics with multiplex and singleplex analysis. *Allergo J Int*. 2015;24:46-53.
11. Williams P et al. Evaluation of a novel automated allergy microarray platform compared with three other allergy test methods. *Clin Exp Immunol*. 2016;184:1-10.
12. Fedenko E et al. Microarray-based IgE serology improves management of severe atopic dermatitis in two children. *Pediatric Allergy and Immunology* 2016;27;645-659.
13. Mari A et al. The IgE-microarray testing in atopic dermatitis: a suitable modern tool for the immunological and clinical phenotyping of the disease. *Curr Opin Allergy Clin Immunol* 11:438-444.
14. Ott H et al. Allergen microarrays: a novel tool for high-resolution IgE profiling in adults with atopic dermatitis *Eur J Dermatol* 2010; 20 (1): 54-61.
15. Choi JS et al. Clinical availability of component-resolved diagnosis using microarray technology in atopic dermatitis. *Ann Dermatol*. 2014;26:437-46.
16. Heaps A et al. The utility of the ISAC allergen array in the investigation of idiopathic anaphylaxis. *Clin Exp Immunol*. 2014 Aug;177(2):483-90.
17. Sastre J et al. How molecular diagnosis can change allergen-specific immunotherapy prescription in a complex pollen area. *Allergy*. 2012 May;67(5):709-11.
18. Asero R et al. Component-resolved diagnosis-assisted prescription of allergen specific immunotherapy: a practical guide. *Eur Ann Allergy Clin Immunol* 2012;44:183-7.
19. Melioli G et al. Potential of molecular based diagnostics and its impact on allergen immunotherapy. *Asthma Research and Practice* 2016 2:9.
20. Kukkonen A K et al. Ara h 2 and Ara h 6 are the best predictors of severe peanut allergy: a double-blind placebo-controlled study. *Allergy* 2015;70:1239-1245.
21. D'Urbano L E et al. Performance of a component-based allergen-microarray in the diagnosis of cow's milk and hen's egg allergy. *Clinical & Experimental Allergy*. 2010;40;1561-1570.
22. Scala E et al. Lipid transfer protein sensitization: reactivity profiles and clinical risk assessment in an Italian cohort. *Allergy* 2015;70:933-943.
23. Nordlund B et al. IgE antibodies to animal-derived lipocalin, kallikrein and secretoglobulin are markers of bronchial inflammation in severe childhood asthma. *Allergy* 2012;67:661-9.
24. Konradsen JR et al. Allergy to Furry animals: New Insights, diagnostic approaches, and challenges. *J Allergy Clin Immunol* 2015;135:616-25.

# Empfohlene Informationsquellen

**allergyai.com** – Website  
der ImmunoDiagnostics von  
Thermo Fisher Scientific

**allergen.org** – Allergenkomponenten-  
datenbank: Systematic allergen nomenclature  
approved by the World Health Organization  
and International Union of Immunological  
Societies (WHO/IUIS) Allergen Nomenclature  
Sub-committee.

Matricardi PM et al. EAACI Molecular  
Allergology User's Guide. Pediatric allergy  
and immunology: official publication of the  
European Society of Pediatric Allergy and  
Immunology. 2016;27 Suppl 23:1-250.

Kleine-Tebbe J and Jakob T Editors:  
Molecular Allergy Diagnostics. Innovation  
for a Better Patient Management. Springer  
International Publishing Switzerland 2017.  
ISBN 978-3-319-42498-9 ISBN 978-3-319-  
42499-6 (eBook), DOI 10.1007/978-3-319-  
42499-6.

Canonica GW et al. A WAO - ARIA - GA<sup>2</sup>LEN  
consensus document on molecular-  
based allergy diagnostics. World Allergy  
Organ J 2013;6(1):17.

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### Literatur:

1. Matricardi PM et al. EAACI Molecular Allergy User's Guide. Pediatric allergy and immunology: official publication of the European Society of Pediatric Allergy and Immunology. 2016;27 Suppl 23:1-250.
2. Kleine-Tebbe J and Jakob T Editors: Molecular Allergy Diagnostics. Innovation for a Better Patient Management. Springer International Publishing Switzerland 2017. ISBN 978-3-319-42498-9 ISBN 978-3-319-42499-6 (eBook), DOI 10.1007/978-3-319-42499-6.
3. Canonica GW et al. A WAO - ARIA - GA<sup>2</sup>LEN consensus document on molecular-based allergy diagnostics. World Allergy Organ J. 2013 Oct 3;6(1):17.
4. Wickman M. When allergies complicate allergies. Allergy 2005;60(S79):14-18.
5. Van Hage M et al. ImmunoCAP assays: Pros and cons in allergology. J Allergy Clin Immunol 2017;140:974-7.

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ImmunoCAP Allergen f13, Peanut; ImmunoCAP Allergen f422, Allergen component rAra h 1 Peanut; ImmunoCAP Allergen f423, Allergen component rAra h 2 Peanut; ImmunoCAP Allergen f424, Allergen component rAra h 3 Peanut; ImmunoCAP Allergen f447, Allergen component rAra h 6 Peanut; ImmunoCAP Allergen f352, Allergen component rAra h 8 Peanut; ImmunoCAP Allergen f427, Allergen component rAra h 9 Peanut; ImmunoCAP Allergen f14, Soybean, ImmunoCAP Allergen f353, Allergen component rGly m 4 PR-10, Soy, ImmunoCAP Allergen f431, Allergen component nGly m 5 Beta-conglycinin, Soy, ImmunoCAP Allergen f432, Allergen component nGly m 6 Glycinin, Soy; ImmunoCAP Allergen f17, Hazel nut; ImmunoCAP Allergen f422, Allergen component rCor a 1 PR-10 Hazel nut; ImmunoCAP Allergen f425, Allergen component rCor a 8, Hazel nut; ImmunoCAP Allergen f440, Allergen component nCor a 9, Hazelnut; ImmunoCAP Allergen f439, Allergen component rCor a 14, Hazelnut; ImmunoCAP Allergen f256, Walnut; ImmunoCAP Allergen f441, Allergen component rJug r 1, Walnut; ImmunoCAP Allergen f442, Allergen component rJug r 3 LTP, Walnut; ImmunoCAP Allergen f202, Cashew nut; ImmunoCAP Allergen f443, Allergen component rAna o 3, Cashew nut; ImmunoCAP Allergen f18, Brazil nut; ImmunoCAP Allergen f354, Allergen component rBer e 1, Brazil nut; ImmunoCAP Allergen f449, Allergen Component rSes i 1, Sesame seed; ImmunoCAP Allergen f49, Apple; ImmunoCAP Allergen f237, Apricot; ImmunoCAP Allergen f95, Peach; ImmunoCAP Allergen f94, Pear; ImmunoCAP Allergen f255, Plum; ImmunoCAP Allergen f20, Almond; ImmunoCAP Allergen f343, Raspberry; ImmunoCAP Allergen f44, Strawberry; ImmunoCAP Allergen f419, Allergen component rPru p 1 PR-10, Peach; ImmunoCAP Allergen f420, Allergen component rPru p 3 LTP, Peach; ImmunoCAP Allergen f421, Allergen component rPru p 4 Profilin, Peach; ImmunoCAP Allergen f454, Allergen Component rPru p 7, Peach; ImmunoCAP Allergen f434, Allergen component rMal d 1 PR-10, Apple; ImmunoCAP Allergen f435, Allergen component rMal d 3 LTP, Apple; ImmunoCAP Allergen f4, Wheat; ImmunoCAP Allergen f98, Gliadin; ImmunoCAP Allergen f416, Allergen component rTri a 19 Omega-5 Gliadin, Wheat; ImmunoCAP Allergen f433, Allergen component rTri a 14 LTP, Wheat; ImmunoCAP Allergen f1, Egg white; ImmunoCAP Allergen f75, Egg yolk; ImmunoCAP Allergen f233, Allergen component nGal d 1 Ovomucoid, Egg; ImmunoCAP Allergen f232, Allergen component nGal d 2 Ovalbumin, Egg; ImmunoCAP Allergen f323, Allergen component nGal d 3 Conalbumin, Egg; ImmunoCAP Allergen k208, Allergen component

nGal d 4 Lysozyme, Egg; ImmunoCAP Allergen f2, Milk; ImmunoCAP Allergen f76, Allergen component nBos d 4 Alpha-lactalbumin, Milk; ImmunoCAP Allergen f77, Allergen component nBos d 5 Beta-lactoglobulin, Milk; ImmunoCAP Allergen e204, Allergen component nBos d 6 BSA, Cow; ImmunoCAP Allergen f78, Allergen component nBos d 8 Casein, Milk; ImmunoCAP Allergen f27, Beef; ImmunoCAP Allergen f26, Pork; ImmunoCAP Allergen f88, Mutton; ImmunoCAP Allergen c74, Gelatin bovine; ImmunoCAP Allergen o215, Component nGal-alpha-1,3-Gal (alpha-Gal) Thyroglobulin, bovine; ImmunoCAP Allergen f24, Shrimp; ImmunoCAP Allergen f23, Crab; ImmunoCAP Allergen f37, Blue mussel; ImmunoCAP Allergen f351, Allergen component rPen a 1 Tropomyosin, Shrimp; ImmunoCAP Allergen d205, Allergen component rDer p 10 Tropomyosin, House dust mite; ImmunoCAP Allergen f3, Fish (cod); ImmunoCAP Allergen f42, Haddock; ImmunoCAP Allergen f41, Salmon; ImmunoCAP Allergen f206, Mackerel; ImmunoCAP Allergen f426, Allergen component rGad c1 Cod; ImmunoCAP Allergen f355, Allergen component rCyp c1 Carp; ImmunoCAP Allergen e1, Cat dander, ImmunoCAP Allergen e94, Allergen component rFel d 1 Cat, ImmunoCAP Allergen e220, Allergen component rFel d 2 Cat serum albumin, ImmunoCAP Allergen e228, Allergen component rFel d 4, Cat, ImmunoCAP Allergen e231, Allergen component rFel d 7, Cat; ImmunoCAP Allergen e5, Dog dander, ImmunoCAP Allergen e101, Allergen component rCan f 1 Dog, ImmunoCAP Allergen e102, Allergen component rCan f 2 Dog, ImmunoCAP Allergen e221, Allergen component nCan f 3 Dog serum albumin, ImmunoCAP Allergen e229, Allergen component rCan f 4, Dog, ImmunoCAP Allergen e226, Allergen component rCan f 5, Dog, ImmunoCAP Allergen e230, Allergen component rCan f 6, Dog; ImmunoCAP Allergen e3, Horse dander; ImmunoCAP Allergen e227, Allergen component rEqu c 1, Horse; ImmunoCAP Allergen d1, House dust mite; ImmunoCAP Allergen d2, House dust mite; ImmunoCAP Allergen d202, Allergen component nDer p 1, House dust mite; ImmunoCAP Allergen d203, Allergen component rDer p 2, House dust mite; ImmunoCAP Allergen d205, Allergen component rDer p 10 Tropomyosin, House dust mite; ImmunoCAP Allergen d209, Allergen component rDer p 23, House dust mite; ImmunoCAP Allergen g2, Bermuda grass; ImmunoCAP Allergen g6, Timothy; ImmunoCAP Allergen g216, Allergen component nCyn d 1 Bermuda grass; ImmunoCAP Allergen g205, Allergen component rPhl p 1 Timothy; ImmunoCAP Allergen g206, Allergen component rPhl p 2 Timothy; ImmunoCAP Allergen g208, Allergen component nPhl p 4 Timothy; ImmunoCAP Allergen g215, Allergen

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component rPhl p 5b Timothy; ImmunoCAP Allergen g209, Allergen component rPhl p 6 Timothy; ImmunoCAP Allergen g210, Allergen component rPhl p 7 Timothy; ImmunoCAP Allergen g211, Allergen component rPhl p 11 Timothy; ImmunoCAP Allergen g212, Allergen component rPhl p 12 Profilin, Timothy; ImmunoCAP Allergen g213, Allergen component rPhl p 1, rPhl p 5b Timothy; ImmunoCAP Allergen g214, Allergen component rPhl p 7, rPhl p 12 Timothy; ImmunoCAP Allergen o214, Allergen component MUXF3 CCD, Bromelain; ImmunoCAP Allergen t3, Common silver birch; ImmunoCAP Allergen t215, Allergen component rBet v 1 PR-10, Birch; ImmunoCAP Allergen t216, Allergen component rBet v 2 Profilin, Birch; ImmunoCAP Allergen t220, Allergen component rBet v 4 Birch; ImmunoCAP Allergen t225, Allergen component rBet v 6 Birch; ImmunoCAP Allergen t221, Allergen component rBet v 2, rBet v 4 Birch; ImmunoCAP Allergen t23, Italian/Mediterranean/Funeral cypress; ImmunoCAP Allergen t222, Arizona cypress; ImmunoCAP Allergen t9, Olive; ImmunoCAP Allergen t11, Maple leaf sycamore, London plane; ImmunoCAP Allergen t226, Allergen component nCup a 1 Cypress; ImmunoCAP Allergen t224, Allergen Component rOle e 1, Olive; ImmunoCAP Allergen t227, Allergen component nOle e 7 LTP, Olive; ImmunoCAP Allergen t240, Allergen Component rOle e 9, Olive; ImmunoCAP Allergen t241, Allergen component rPla a 1, Maple leaf sycamore, London plane; ImmunoCAP Allergen w6, Mugwort; ImmunoCAP Allergen w21, Wall pellitory; ImmunoCAP Allergen w9, Plantain (English), Ribwort; ImmunoCAP Allergen w11, Saltwort (prickly), Russian thistle; ImmunoCAP Allergen w230, Allergen component nAmb a 1 Ragweed; ImmunoCAP Allergen w231, Allergen component nArt v 1 Mugwort; ImmunoCAP Allergen w233, Allergen component nArt v 3 LTP, Mugwort; ImmunoCAP Allergen w211, Allergen component rPar j 2 LTP, Wall pellitory; ImmunoCAP Allergen w234, Allergen component rPla l 1, Plantain; ImmunoCAP Allergen w232, Allergen component nSal k 1 Saltwort; ImmunoCAP Allergen m6, Alternaria alternate; ImmunoCAP Allergen m229, Allergen component rAlt a 1, Alternaria alternata; ImmunoCAP Allergen m3, Aspergillus fumigatus; ImmunoCAP Allergen m218, Allergen component rAsp f 1 Aspergillus fumigatus; ImmunoCAP Allergen m219, Allergen component rAsp f 2 Aspergillus fumigatus; ImmunoCAP Allergen m220, Allergen component rAsp f 3 Aspergillus fumigatus; ImmunoCAP Allergen m221, Allergen component rAsp f 4 Aspergillus fumigatus; ImmunoCAP Allergen m222, Allergen component rAsp f 6 Aspergillus fumigatus; ImmunoCAP Allergen i1, Honey bee venom; ImmunoCAP Allergen i3, Common wasp venom (Yellow

jacket); ImmunoCAP Allergen i4, Paper wasp venom; ImmunoCAP Allergen i208, Allergen component rApi m 1 Phospholipase A2, Honey bee; ImmunoCAP Allergen i214, Allergen component rApi m 2, Honey bee; ImmunoCAP Allergen i215, Allergen component rApi m 3, Honey bee; ImmunoCAP Allergen i216, Allergen component rApi m 5, Honey bee; ImmunoCAP Allergen i217, Allergen component rApi m 10, Honey bee; ImmunoCAP Allergen i211, Allergen component rVes v 1 Phospholipase A1, Common wasp; ImmunoCAP Allergen i209, Allergen component rVes v 5 Common wasp; ImmunoCAP Allergen i210, Allergen component rPol d 5 European Paper wasp; ImmunoCAP Allergen k82, Latex; ImmunoCAP Allergen k218, Allergen component rHev b 5 Latex; ImmunoCAP Rare Allergen k215, Allergen component rHev b 1 Latex; ImmunoCAP Rare Allergen k217, Allergen component rHev b 3 Latex; ImmunoCAP Rare Allergen k220, Allergen component rHev b 6.02 Latex; ImmunoCAP Rare Allergen k221, Allergen component rHev b 8 Profilin, Latex; ImmunoCAP Rare Allergen k224, Allergen component rHev b 11 Latex.

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© 2021 Thermo Fisher Scientific Inc. Alle Rechte vorbehalten. Alle Warenzeichen sind das Eigentum von Thermo Fisher Scientific und seiner Tochtergesellschaften, falls nicht anders angegeben.
Rechtmäßiger Hersteller: Phadia AB, Uppsala, Schweden 75163.AL.EU2.DE.v1.21 84213210
**Thermo Fisher Diagnostics GmbH**, Munzinger Str. 7, D-79111 Freiburg, Tel. +49 761 47 805 0, Fax +49 761 47 805 338
**Thermo Fisher Diagnostics Austria GmbH**, Dresdner Str. 89, A-1200 Wien, Tel. +43 1 270 20 20, Fax +43 1 270 20 20 20
**Thermo Fisher Diagnostics AG**, Senneweidstr. 46, CH-6312 Steinhausen, Tel. +41 43 343 40 50, Fax +41 43 343 40 51