Rheumatoide Arthritis

Gut zu Wissen

Rheumatoide Arthritis (RA)

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Referenzen

- Aletaha D, Blum S. Therapeutic implications of autoantibodies in rheumatoid arthritis. *RMD Open*. 2016;2(1):e000009.
- Littlejohn EA, Monrad SU. Early Diagnosis and Treatment of Rheumatoid Arthritis. *Prim Care*. 2018;45(2):237-55.
- de Brito Rocha S, Baldo DC, Andrade LEC. Clinical and pathophysiologic relevance of autoantibodies in rheumatoid arthritis. *Adv Rheumatol*. 2019;59(1):2.
- Smolen JS, Landewe R, Bijlsma J, Burmester G, Chatzidionysiou K, Dougados M, et al. EULAR recommendations for the management of rheumatoid arthritis with synthetic and biological disease-modifying antirheumatic drugs: 2016 update. *Ann Rheum Dis*. 2017;76(6):960-77.
- Bizzaro N. Antibodies to citrullinated peptides: a significant step forward in the early diagnosis of rheumatoid arthritis. *Clin Chem Lab Med*. 2007;45(2):150-7.
- Aletaha D, Neogi T, Silman AJ, Funovits J, Felson DT, Bingham CO, 3rd, et al. 2010 Rheumatoid arthritis classification criteria: an American College of Rheumatology/European League Against Rheumatism collaborative initiative. *Arthritis Rheum*. 2010;62(9):2569-81.
- Santos-Moreno P, Villarreal-Peralta L, Gomez-Mora D, Castillo E, Malpica A, Castro C, et al. Osteoarthritis is the Most Frequent Cause of Rheumatoid Arthritis Misdiagnosis in a Rheumatology Center. *Osteoarthritis and Cartilage*. 2017;25:S219.
- van Venrooij WJ, van Beers JJ, Pruijn GJ. Anti-CCP antibodies: the past, the present and the future. *Nat Rev Rheumatol*. 2011;7(7):391-8.
- Mathsson Alm L, Fountain DL, Cadwell KK, Madrigal AM, Gallo G, Poorafshar M. The performance of anti-cyclic citrullinated peptide assays in diagnosing rheumatoid arthritis: a systematic review and meta-analysis. *Clin Exp Rheumatol*. 2018;36(1):144-52.
- Simundic AM. Measures of Diagnostic Accuracy: Basic Definitions. *EJIFCC*. 2009;19(4):203-11.
- Barhamain AS, Magliah RF, Shaheen MH, Munassar SF, Falemban AM, Alshareef MM, et al. The journey of rheumatoid arthritis patients: a review of reported lag times from the onset of symptoms. *Open Access Rheumatol*. 2017;9:139-50.
- Johnson PM, Faulk WP. Rheumatoid factor: its nature, specificity, and production in rheumatoid arthritis. *Clin Immunol Immunopathol*. 1976;6(3):414-30.
- Williams DG, Moyes SP, Mageed RA. Rheumatoid factor isotype switch and somatic mutation variants within rheumatoid arthritis synovium. *Immunology*. 1999;98(1):123-36.
- Falkenburg WJJ, van Schaardenburg D. Evolution of autoantibody responses in individuals at risk of rheumatoid arthritis. *Best Pract Res Clin Rheumatol*. 2017;31(1):42-52.
- Rantapaa-Dahlqvist S, de Jong BA, Berglin E, Hallmans G, Wadell G, Stenlund H, et al. Antibodies against cyclic citrullinated peptide and IgA rheumatoid factor predict the development of rheumatoid arthritis. *Arthritis Rheum*. 2003;48(10):2741-9.
- Jonsson T, Thorsteinsson J, Valdimarsson H. Elevation of only one rheumatoid factor isotype is not associated with increased prevalence of rheumatoid arthritis—a population based study. *Scand J Rheumatol*. 2000;29(3):190-1.
- Halldorsdottir HD, Jonsson T, Thorsteinsson J, Valdimarsson H. A prospective study on the incidence of rheumatoid arthritis among people with persistent increase of rheumatoid factor. *Ann Rheum Dis*. 2000;59(2):149-51.
- Jaskowski TD, Hill HR, Russo KL, Lakos G, Szekanecz Z, Teodorescu M. Relationship between rheumatoid factor isotypes and IgG anti-cyclic citrullinated peptide antibodies. *J Rheumatol*. 2010;37(8):1582-8.
- Ingegnoli F, Castelli R, Gualtierotti R. Rheumatoid factors: clinical applications. *Dis Markers*. 2013;35(6):727-34.
- Jonsson T, Steinsson K, Jonsson H, Geirsson AJ, Thorsteinsson J, Valdimarsson H. Combined elevation of IgM and IgA rheumatoid factor has high diagnostic specificity for rheumatoid arthritis. *Rheumatol Int*. 1998;18(3):119-22.
- Hamilton RG. Methods (In Vitro and In Vivo): Nephelometry and Turbidimetry. In: Mackay IR, Rose NR, Ledford DK, Lockey RF, editors. *Encyclopedia of Medical Immunology: Allergic Diseases*. New York, NY: Springer New York; 2014. p. 484-6.
- Bas S, Perneger TV, Kunzle E, Vischer TL. Comparative study of different enzyme immunoassays for measurement of IgM and IgA rheumatoid factors. *Ann Rheum Dis*. 2002;61(6):505-10.
- de Angelis V, Meroni PL. *Autoantibodies*. 2nd ed. Amsterdam ; Boston: Elsevier; 2007. xxiv, 838 p.
- Mikulic TR, O'Dell JR, Stoner JA, Parrish LA, Arend WP, Norris JM, et al. Association of rheumatoid arthritis treatment response and disease duration with declines in serum levels of IgM rheumatoid factor and anti-cyclic citrullinated peptide antibody. *Arthritis Rheum*. 2004;50(12):3776-82.
- Bobbio-Pallavicini F, Caporali R, Alpini C, Avalle S, Epis OM, Klersy C, et al. High IgA rheumatoid factor levels are associated with poor clinical response to tumour necrosis factor alpha inhibitors in rheumatoid arthritis. *Ann Rheum Dis*. 2007;66(3):302-7.
- Sieghart D, Platzer A, Studenic P, Alasti F, Grundhuber M, Swiniarski S, et al. Determination of Autoantibody Isotypes Increases the Sensitivity of Serodiagnostics in Rheumatoid Arthritis. *Front Immunol*. 2018;9:876.

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© 2020 Thermo Fisher Scientific Inc. Alle Rechte vorbehalten. Alle Warenzeichen sind das Eigentum von Thermo Fisher Scientific und seiner Tochtergesellschaften, falls nicht anders angegeben. Rechtmäßiger Hersteller: Phadia AB, Uppsala, Schweden.
Thermo Fisher Diagnostics GmbH, Munzinger Str. 7, D-79111 Freiburg, Tel. +49 761 47 805 0, Fax +49 761 47 805 338
Thermo Fisher Diagnostics Austria GmbH, Dresdner Str. 89, A-1200 Wien, Tel. +43 1 2702020, Fax +43 1 270202020
Thermo Fisher Diagnostics AG, Sennweidstr. 46, CH-6312 Steinhausen, Tel. +41 43 343 4050, Fax +41 43 343 4051