




Setting the standard

ImmunoCAP™ Specific IgE Stinging Insect Allergen Components*

Whole allergen and allergen component testing provides additional information that can help diagnose allergy more precisely, allowing for a more comprehensive management plan.¹

Testing with stinging insect components can help to:¹⁻⁹

-  Discriminate between species specific sensitization and cross-reactivity⁴
-  Identify culprit venom(s)⁵
-  Facilitate accurate prescription of venom immunotherapy (VIT)⁵

Up to
50%

of venom allergic patients test positive for both bee and wasp venom.²

Component resolved IgE tests using recombinant venom allergens may improve specificity³—increasing the likelihood of successful venom immunotherapy.⁴

Stinging Insect Profile

Hymenoptera venom allergy profiles may consist of five whole allergens, eight components, Tryptase and CCD.⁶




i1. Honeybee

i3. Common Wasp (Yellow Jacket)

i4. Paper Wasp

i2. White Faced Hornet†

i5. Yellow Hornet†


	Honeybee Apis mellifera 	Yellow Jacket Vespula vulgaris 	Paper Wasp Polistes dominula 
Differentiating marker allergen	Api m 1 Phospholipase A2 Api m 3 Acid phosphatase Api m 10 Icarapin	<div style="border: 1px dashed blue; padding: 5px;"> Ves v 1 Phosphlipase A1 Ves v 5 Antigen 5 </div>	
	Api m 2 Hyaluronidase Api m 5 Dipeptidyl peptidase		

Adapted from S. Blank, M. B. Bilo, M. Ollert Component-resolved diagnostics to direct in venom immunotherapy: Important steps towards precision medicine Clin Exp Allergy. 2018;48:p357. 2018

†Components not currently available

Management Considerations⁶⁻⁸

Results should be interpreted in the context of a patient's clinical symptoms and history.
CCD-Bromelain (MUXF3) should be measured to assess cross-reactivity between species.⁹



“**Tryptase should be analyzed in patients with a history of a severe sting reaction and before starting venom SIT.**”^{3,10-14}

EAACI, AAAAI, WAO, ICON

ImmunoCAP™ Tryptase is FDA cleared for in vitro diagnostic use as an aid in the clinical diagnosis of patients with a suspicion of systemic mastocytosis in conjunction with other clinical findings.

*Official product names of allergen components mentioned within this document: ImmunoCAP Allergen i1, Honey bee venom; ImmunoCAP Allergen i208, Allergen Component rApi m 1; ImmunoCAP Allergen i214, Allergen Component rApi m 2; ImmunoCAP Allergen i215, Allergen Component rApi m 3; ImmunoCAP Allergen i216, Allergen Component rApi m 5; ImmunoCAP Allergen i217, Allergen Component rApi m 10; ImmunoCAP Allergen i4, Paper wasp; ImmunoCAP Allergen i210, Allergen Component rPol d 5, European Paper wasp; ImmunoCAP Allergen i3, Common wasp venom (Yellow jacket); ImmunoCAP Allergen i211, Allergen Component rVes v 1, Common wasp; ImmunoCAP Allergen i209, Allergen Component rVes v 5, Common wasp; ImmunoCAP Allergen i5, Yellow hornet venom; ImmunoCAP Allergen i2, White-faced hornet venom

References

- Kleine-Tebbe J, Jappe U. Molecular allergy diagnostic tests: development and relevance in clinical practice. *Allergologie select.* 2017;1 (2):169-1893.
- Spillner E. et al., Hymenoptera allergens: from venom to "venome". *Frontiers in immunology* 2014; 5: 1-7.
- Golden, et al Stinging insect hypersensitivity A practice parameter update 2016 *Ann Allergy Asthma Immunol* 118 (2017) p28-54
- Matricardi PM, et al EAACI Molecular Allergy User's Guide. *Pediatr Allergy Immunol* 2016; 27: (suppl23): p157. 2016
- Silke C, Hofmann, MD, et al Added value of IgE detection to rApi m 1 and rVes v 5 in patients with Hymenoptera venom allergy *J ALLERGY CLIN IMMUNOL VOLUME 127, NUMBER 1 265-267 2011*
- Ansotegui I. A WAO — ARIA — GA2LEN consensus document on molecular-based allergy diagnosis (PAMD®): Update 2020. *World Allergy Organization Journal.* 2020;13(2):100091.
- S. Blank, M. B. Bilo, M. Ollert Component-resolved diagnostics to direct in venom immunotherapy: Important steps towards precision medicine *Clin Exp Allergy.* 2018;48:p357. 2018
- Barber D, Diaz-Perales A, Escribese M, Kleine-Tebbe J, Matricardi P, Ollert M, Santos A, Sastre J. Molecular allergology and its impact in specific allergy diagnosis and therapy. *Allergy.* 2021;76(12):3642-3658.
- Irene Mittermann, PhD, Mihaela Zidarn, MD, Mira Silar, P7hD, Zora Markovic-Housley, PhD, Werner Aberer, MD, Peter Korosec, PhD, Mitja Kosnik, PhD, MD, and Rudolf Valenta, MD. Recombinant allergen- based IgE testing to distinguish bee and wasp allergy. *J ALLERGY CLIN IMMUNOL VOLUME 125, NUMBER 6 p1300-1307. 2010*
- Bonifazi F. et al & EAACI Interest Group on Insect Venom Hypersensitivity, Prevention and treatment of hymenoptera venom allergy: guidelines for clinical practice. *Allergy* 2005; 60: 1459-1470.
- Bil B. et al & EAACI Interest Group on Insect Venom Hypersensitivity., Diagnosis of Hymenoptera venom allergy. *Allergy* 2005; 60: 1339-49
- Simons FE. et al., International concensus on (ICON) anaphylaxis. *World Allergy Organ J.* 2014 May 30;7(1):9.
- Simons FE. et al., World Allergy Organization Anaphylaxis Guidelines:2013 update of the evidence base. *Int Arch Allergy Immunol.*2013;162(3):193-204.
- Cox L. et al., Allergen Immunotherapy: A practice parameter third update. *J Allergy Clin Immunol* 2011;127(1):s1-s55.

Learn more at thermofisher.com/allergencomponents