# ReCAPtulando

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#### List of full products names - Main Allergens

ImmunoCAP™ allergen products: ImmunoCAP allergen f96, Avocado; ImmunoCAP allergen f210, Pineapple; ImmunoCAP allergen f225, Pumpkin; ImmunoCAP allergen f47, Garlic; ImmunoCAP allergen f255, Plum; ImmunoCAP allergen f20, Almond; ImmunoCAP allergen f13, Peanut; ImmunoCAP allergen f9, Rice; ImmunoCAP allergen f40, Tuna; ImmunoCAP allergen f7, Oat; ImmunoCAP allergen f17, Hazel nut; ImmunoCAP allergen f92, Banana; ImmunoCAP allergen f35, Potato; ImmunoCAP allergen f300, Goat milk; ImmunoCAP allergen f93, Cacao; ImmunoCAP allergen f24, Shrimp; ImmunoCAP allergen f23, Crab; ImmunoCAP allergen f202, Cashew nut; ImmunoCAP allergen f18, Brazil nut; ImmunoCAP allergen f48, Onion; ImmunoCAP allergen f31, Carrot; ImmunoCAP allergen f242, Cherry; ImmunoCAP allergen f6, Barley; ImmunoCAP allergen f36, Coconut; ImmunoCAP Rare allergen f340, Cochineal extract (Carmine red); ImmunoCAP allergen f12, Pea; ImmunoCAP allergen f214, Spinach; ImmunoCAP allergen f15, White bean; ImmunoCAP allergen f83, Chicken; ImmunoCAP allergen f10, Sesame seed; ImmunoCAP allergen f79, Gluten; ImmunoCAP allergen f84, Kiwi fruit; ImmunoCAP allergen f80, Lobster; ImmunoCAP allergen f33, Orange; ImmunoCAP allergen f2, Milk; ImmunoCAP allergen f208, Lemon; ImmunoCAP allergen f58, Pacific squid; ImmunoCAP allergen f49, Apple; ImmunoCAP allergen f293, Papaya; ImmunoCAP allergen f91, Mango; ImmunoCAP allergen f247, Honey; ImmunoCAP allergen f87, Melon; ImmunoCAP allergen f37, Blue mussel; ImmunoCAP allergen f8, Maize, Corn; ImmunoCAP allergen f44, Strawberry; ImmunoCAP allergen f256, Walnut; ImmunoCAP allergen f245, Egg; ImmunoCAP allergen f1, Egg white; ImmunoCAP allergen f75, Egg yolk; ImmunoCAP allergen f3, Fish (cod); ImmunoCAP Allergen f94, Pear; ImmunoCAP Allergen f95, Peach; ImmunoCAP allergen f280, Black pepper; ImmunoCAP allergen f279, Chili pepper; ImmunoCAP allergen f263, Green pepper (unripe seed); ImmunoCAP allergen f218, Paprika, Sweet pepper; ImmunoCAP allergen f26, Pork; ImmunoCAP allergen f81, Cheese, cheddar type; ImmunoCAP allergen f82, Cheese, mold type; ImmunoCAP allergen f14, Soybean; ImmunoCAP allergen f41, Salmon; ImmunoCAP allergen f61, Sardine, Pilchard; ImmunoCAP allergen f259, Grape; ImmunoCAP allergen f259, Grape; ImmunoCAP allergen f27, Beef; ImmunoCAP allergen c8, Chlorhexidine; ImmunoCAP allergen c261, Pholcodine; ImmunoCAP allergen c71, Insulin bovine; ImmunoCAP allergen c70, Insulin porcine; ImmunoCAP allergen c73, Insulin human; ImmunoCAP allergen c260, Morphine; ImmunoCAP allergen c1, Penicilloyl G; ImmunoCAP allergen c2, Penicilloyl V; ImmunoCAP allergen m6, *Alternaria alternata*; ImmunoCAP allergen m3, *Aspergillus fumigatus*; ImmunoCAP allergen m5, *Candida albicans*; ImmunoCAP allergen m2, *Cladosporium herbarum*; ImmunoCAP allergen m209, *Penicillium glabrum*; ImmunoCAP allergen m80, *Staphylococcal enterotoxin A*; ImmunoCAP allergen m81, *Staphylococcal enterotoxin B*; ImmunoCAP allergen m223, *Staphylococcal enterotoxin C*; ImmunoCAP allergen i206, Cockroach, American; ImmunoCAP allergen i6, Cockroach, German; ImmunoCAP allergen i70, Fire ant; ImmunoCAP allergen i204, Horse fly; ImmunoCAP allergen i71, Mosquito; ImmunoCAP allergen i1, Honey bee venom; ImmunoCAP allergen i4, Paper wasp venom; ImmunoCAP allergen t9, Olive; ImmunoCAP allergen f423, allergen component rAra h 2 Peanut; ImmunoCAP allergen f424, allergen component rAra h 3 Peanut; ImmunoCAP allergen g2, Bermuda grass; ImmunoCAP allergen g5, Rye-grass; ImmunoCAP allergen g6, Timothy; ImmunoCAP allergen g8, Meadow grass, Kentucky blue; ImmunoCAP allergen g10, Johnson grass; ImmunoCAP allergen g17, Bahia grass; ImmunoCAP allergen d201, House dust mite; ImmunoCAP allergen d2, House dust mite; ImmunoCAP allergen d3, House dust mite; ImmunoCAP allergen d1, House dust mite; ImmunoCAP allergen h2, House dust (Hollister-Stier Labs.); ImmunoCAP allergen e5, Dog dander; ImmunoCAP allergen e85, Chicken feathers; ImmunoCAP allergen e1, Cat dander; ImmunoCAP allergen e4, Cow dander; ImmunoCAP allergen o1, Cotton, crude fibers; ImmunoCAP allergen k82, Latex; ImmunoCAP allergen o201, Tobacco leaf; ImmunoCAP allergen d202, allergen component rDer p 1, House dust mite; ImmunoCAP allergen d203, allergen component rDer p 2, House dust mite; ImmunoCAP allergen d205, allergen component rDer p 10 Tropolymyosin, House dust mite; ImmunoCAP allergen f76, allergen component nBos d 4 Alpha-lactalbumin, Milk; ImmunoCAP allergen f77, allergen component nBos d 5 Beta-lactoglobulin, Milk; ImmunoCAP allergen f78, allergen component nBos d 8 Casein, Milk; ImmunoCAP allergen f323, allergen component nGal d 3 Conalbumin, Egg; ImmunoCAP Rare allergen f333, Linseed; ImmunoCAP allergen f98, Gliadin; ImmunoCAP allergen f416, allergen component rTri a 19 Omega-5 Gliadin, Wheat; ImmunoCAP allergen f433, allergen component rTri a 14 LTP, Wheat; ImmunoCAP allergen f431, allergen component nGly m 5 Beta-conglycinin, Soy; ImmunoCAP allergen f432, allergen component nGly m 6 Glycinin, Soy; ImmunoCAP allergen f353, allergen component rGly m 4 PR-10, Soy; ImmunoCAP allergen f422, allergen component rAra h 1 Peanut; ImmunoCAP allergen f423, allergen component rAra h 2 Peanut; ImmunoCAP allergen f424, allergen component rAra h 3 Peanut; ImmunoCAP allergen f352, allergen component rAra h 8 PR-10, Peanut; ImmunoCAP allergen f443, allergen component rAna o 3, Cashew nut; ImmunoCAP allergen f354, allergen component rBer e 1, Brazil nut; ImmunoCAP allergen f428, allergen component rCor a 1 PR-10 Hazel nut; ImmunoCAP allergen f439, allergen component rCor a 14, Hazelnut; ImmunoCAP allergen f25, Tomato; ImmunoCAP allergen f440, allergen component nCor a 9, Hazelnut; ImmunoCAP allergen f441, allergen component rJug r 1, Walnut; ImmunoCAP allergen f442, allergen component rJug r 3 LTP, Walnut; ImmunoCAP allergen f355, allergen component rCyp cl Carp; ImmunoCAP allergen f426, Allergen component rGad cl Cod; ImmunoCAP allergen f351, allergen component rPen a 1 Tropolymyosin, Shrimp; ImmunoCAP allergen f419, allergen component rPru p 1 PR-10, Peach; ImmunoCAP allergen f420, allergen component rPru p 3 LTP, Peach; ImmunoCAP allergen f421, allergen component rPru p 4 Profilin, Peach; ImmunoCAP allergen e204, allergen component nBos d 6 BSA, Cow; ImmunoCAP rare allergen e222, allergen component nSus s PSA, Swine; ImmunoCAP allergen e101, allergen component rCan f 1 Dog; ImmunoCAP allergen e102, allergen component rCan f 2 Dog; ImmunoCAP allergen e221, allergen component nCan f 3 Dog serum albumin; ImmunoCAP allergen e226, allergen component rCan f 5, Dog; ImmunoCAP allergen e94, allergen component rFel d 1 Cat; ImmunoCAP allergen e220, allergen component rFel d 2 Cat serum albumin; ImmunoCAP allergen e228, allergen component rFel d 4, Cat; ImmunoCAP allergen e227, allergen component rEqu c 1, Horse; ImmunoCAP rare allergen k215, allergen component rHev b 1 Latex; ImmunoCAP rare allergen k217, allergen component rHev b 3 Latex; ImmunoCAP allergen k218, allergen component rHev b 5 Latex; ImmunoCAP rare allergen k220, allergen component rHev b 6.02 Latex; ImmunoCAP rare allergen k221, allergen component rHev b 8 Profilin, Latex; ImmunoCAP rare allergen k224, allergen component rHev b 11 Latex; ImmunoCAP allergen i208, allergen component rApi m 1 Phospholipase A2, Honey bee; ImmunoCAP allergen i210, allergen component rPol d 5 European Paper wasp; ImmunoCAP allergen i211, allergen component rVes v 1 Phospholipase A1, Common wasp; ImmunoCAP allergen i209, allergen component rVes v 5 Common wasp; ImmunoCAP Allergen m218, Allergen component rAsp f 1 *Aspergillus fumigatus*; ImmunoCAP allergen m229, allergen component rAlt a 1, *Alternaria alternata*; ImmunoCAP allergen o215, Component nGal-alpha-1,3-Gal (alpha-Gal) Thyroglobulin, bovine; ImmunoCAP allergen o214, allergen component MUXF3 CCD, Bromelain.

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
## Referências

1. Thanh D, Dang Rachel L, Peters Jennifer J, Koplin Shyamali C, Dharmage Lyle C, Gurrin Anne Louise, Ponsonby David J, Martino Melanie, Neeland Mimi L, K. Tang Katrina J, Allen. Egg allergen specific IgE diversity predicts resolution of egg allergy in the population cohort HealthNuts. *Allergy*. 2019 Feb;74(2):318-326. doi: 10.1111/all.13572. Epub 2018 Nov 19.
2. Ando H, Movérare R, Kondo Y, et al. Utility of ovomucoid-specific IgE concentrations in predicting symptomatic egg allergy. *J Allergy Clin Immunol*. 2008;122(3):583-588.
3. Benhamou AH, Caubet JC, Eigenmann PA, et al. State of the art and new horizons in the diagnosis and management of egg allergy. *Allergy*. 2010;65(3): 283-289.
4. Shin M, Han Y, Ahn K. The influence of the time and temperature of heat treatment on the allergenicity of egg white proteins. *Allergy Asthma Immunol Res*. 2013;5(2):96-101.
5. Everberg H, Brostedt P, Oman H, Bohman S, Movérare R. Affinity purification of egg-white allergens for improved component-resolved diagnostics. *Int Arch Allergy Immunol* 2011;154:33-41.
6. Clei G, Brusca I, Scala E, Villata D, Pastorello E, Farioli L, et al. Allergy 2019;74(9):1787-9. House dust mite allergy in Italy – Diagnostic and clinical relevance of Der p 23 (and of minor allergens): A real-life multicenter study.
7. Jiménez-Feijoo R, Pascal M, Moya R, Riggioni C, Domínguez O, Lózano J, et al. Molecular diagnosis in house dust mite allergic patients suggests clinical relevance of Der p 23 in asthmatic children. *J Investig Allergol Clin Immunol* 2019; in press.
8. Weghofer M, Resch Y, Casset A, Kneidinger M, Kopec J, et al. Identification of Der p 23, a peritrophin-like protein, as a new major Dermatophagoides pteronyssinus allergen associated with the peritrophic matrix of mite fecal pellets. *J Immunol*. 2013;190(7):3059-67.
9. Simpson A, Soderstrom L, Ahlstedt S, Murray CS, Woodcock A, Custovic A. IgE antibody quantification and the probability of wheeze in preschool children. *The Journal of allergy and clinical immunology*. 2005;116(4):744-9.
10. Soto-Quiros M, Avila L, Platts-Mills TA, Hunt JF, Erdman DD, Carper H, et al. High titers of IgE antibody to dust mite allergen and risk for wheezing among asthmatic children infected with rhinovirus. *The Journal of allergy and clinical immunology*. 2012;129(6):1499- 505 e5.
11. Asero R. *Eur Ann Allergy Clin Immunol*. 2012;44(5):183-7.
12. Schmid-Grendelmeier P. *Hautarzt*. 2010;61(11):946-53.
13. Weghofer M, et al. *J Immunol*. 2013;190(7):3059-67.
14. Posa D, et al. *J Allergy Clin Immunol*. 2017;139:541-94.4
15. Resch Y, et al. *J Allergy Clin Immunol*. 2015;136:1083-91.

Saiba mais em: [allergyinsider.com](http://allergyinsider.com)

**ThermoFisher**  
SCIENTIFIC