



Phadia 250 16-position sample rack

Improve operating efficiency and walk away time

The new 16-position sample rack increases loading capacity of the Phadia™ 250 instrument and extends walk away time without making any changes to your laboratory set up.

Increased sample loading capacity

Load up to 80 sample tubes simultaneously on your Phadia 250 instrument, a 60% increase in capacity

Flexible handling

Load a mix of 10-position and 16-position racks at the same time to suit your workflow needs

Quick and easy set-up

Each rack is uniquely identified by the Phadia™ Prime software



Order information

	Phadia 250 10-position sample rack	Phadia 250 16-position sample rack
Article number	12-4906-80	12-4907-73
Number of racks/box	5 pcs	5 pcs
Sample rack positions	10	16
Tube diameter	10-17 mm	10-13.0 mm

Additional information

- Phadia 250 10-position sample racks are delivered with new Phadia 250 instruments and available for ordering as accessory
- Phadia 250 16-position sample racks are available for ordering as accessory

Find out more at thermofisher.com/phadia