

Order information for ImmunoCAP™ Allergen f449, Allergen Component rSes i 1 Sesame seed

Product description

rSes i 1, 2S Albumin storage protein
Sesame seed (*Sesamum indicum*)

Allergen code

f449

Article number

14-6109-01

Barcode

E7M



Available
from
July 1st, 2021

Note! The test will be included in the instrument software Phadia™ Prime software version 2.6.
To run the test before the update, the article needs to be manually added in the software.

Improved sesame allergy diagnosis using rSes i 1

A new ImmunoCAP™ Allergen Component test providing:



Aid in the diagnosis and management of patients with primary sesame allergy.¹⁻⁵



Help assessing risk for severe events by providing quantitative sIgE levels to Ses i 1.¹⁻⁵



Higher clinical specificity than skin prick test and whole allergen sIgE test with component-resolved diagnostics.²⁻⁴

1. Maruyama N et al. Clin Exp Allergy. 2016;46(1):163-71. 2. Yanagida N et al. J Allergy Clin Immunol Pract. 2019;7(6):2084-86. 3. Saf S et al. J Allergy Clin Immunol Pract. 2020;8(5):1681-88. 4. Goldberg MR et al. Pediatr Allergy Immunol. 2021 May 5. doi: 10.1111/pai.13533. Online ahead of print. 5. Nachshon L et al. J Allergy Clin Immunol Pract. 2019;7:2775-81.

Find out more at thermofisher.com/phadia

© 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. 175251.AL.ALR1.EN.v1.21



ThermoFisher
SCIENTIFIC